# CEDIA™ Theophylline II Assay

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1. Rall TW, The xanthines. In: Gilman, A.G., Goodman, L.S. and Gilman, A. eds: The Pharmacological Basis of Therapeutics. New York: MacMillan Publishing Company, 1980: 592-607.
2. Hendeles L, Weinberger MM: Theophylline therapeutic use and serum concentration monitoring. In: Taylor WJ, Finn AL, eds: Individualizing Drug Therapy: Practical Applications of Drug Monitoring. New York: Gross Townsend Frank, Inc., 1981: vol 1, pp 31-66.
3. Aranda JV, Sitar DS, Parsons WD, Loughnan PM, Neims AH: Pharmacokinetic aspects of theophylline in premature newborns. New England Journal of Medicine 1976; 295: 413-416.
4. Ogilvie RI: Clinical pharmacokinetics of theophylline. Clinical Pharmacokinetics 1978; 3:267-293.
5. Jacobs MH, Senior RM, Kessler G: Clinical experience with theophylline. Relationships between dosage, serum concentration, and toxicity. Journal of the American Medical Association 1976;235: 1983-1986.
6. Weinberger MW, Matthay RA, Ginchansky EJ, Chidsey CA, Petty TL: Intravenous aminophylline dosage. Use of serum theophylline measurement for guidance. Journal of the American Medical Association 1976;235: 2110-2113.
7. Zwillich CW, Sutton FD, Neff TA, Cohn WM, Matthay RA, Weinberger MM: Theophylline-induced seizures in adults. Correlation with serum concentrations. Annals of Internal Medicine 1975;82: 784-787.
8. Piafsky KM, Ogilvie RI: Drug therapy. Dosage of theophylline in bronchial asthma. New England Journal of Medicine 1975; 292:1218-1222.
9. Leung P, Kalisker A, Bell TD: Variation in theophylline clearance rate with time in chronic childhood asthma. Journal of Allergy and Clinical Immunology 1977;59: 440-444.
10. Henderson DR, Friedman SB, Harris JD, Manning WB, Zoccoli MA: CEDIA, a New Homogeneous Immunoassay System. Clin Chem 1986;32(9): 1637-1641.
11. Data on traceability are on file at Microgenics Corporation, a part of Thermo Fisher Scientific.
12. Mitenko, P.A. and Ogilvie, R.I.: Rational intravenous doses of theophylline. New England Journal of Medicine 1973b;289: 600-603.
13. Buelow KB, Larsson H, Leideman T: Plasma theophylline level and ventilatory function in chronic obstructive pulmonary disease during prolonged oral treatment with choline theophyllinate. European Journal of Clinical Pharmacology 1975;8: 119-123.
14. Weinberger MM, Bronsky EA: Evaluation of oral bronchodilator therapy in asthmatic children. Journal of Pediatrics 1974; 84:421-427.
15. Breiner R, McComb Lewis S, Wong, SHY, Marzouk N, Kapke GF: Positive interference with immunoassay of theophylline in serum of uremics, letters to the editor. Clin. Chem. 1985;31:1575-1577.
16. Nelson KM, Mathews SE, Bowers LD: Theophylline concentrations may be falsely high in serum or uremic patients (letter to the editors). Clin Chem 1983;29: 2125-2126.
17. Nicot G, Charnes JP, Lachatre G, Sautereau D, Valette JP, Eichler E, Leroux-Robert C: Theophylline toxicity risks and chronic renal failure. Int J Clin Pharmacol 1989; 27:398-401.
18. Opein KE, Ainaridi V, Raisys VA, Smith CM, Messenfer LJ: Increase in apparent theophylline concentration in the serum of two uremic patients as measured by some immunoassay methods (caused by 1,3-dimethyluric acid?) (letters to the editors). Clin Chem 1983;29:1698-1699.
19. Patel JA, Clayton LT, LeBel CP, McClatchey KD: Abnormal theophylline levels in plasma by fluorescence polarization immunoassay in patients with renal disease. Ther. Drug Monit. 1984;6:458-460.
20. Data on file at Microgenics Corporation, a part of Thermo Fisher Scientific.

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B-R-A-H-M-S GmbH  
Neuendorfstrasse 25  
16761 Hennigsdorf, Germany

若要取得仿單更新，請前往：  
[www.thermofisher.com/diagnostics](http://www.thermofisher.com/diagnostics)

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