# Thermo Scientific Nalgene Carboys

## Autoclave Recommendations for Thermo Scientific™ Nalgene™ Carboys (EN)

- The recommended autoclaving cycle is 121°C, 1 bar and/or 250°F, 15 psig for 15-20 minutes
- Avoid stacking carboys or using autoclave baskets with other objects on top
- Do not secure the openings of the carboy with aluminum foil, blue Steriwrap, gauze, cotton or tape
- Do not autoclave with the closure or closed spigot in place (open spigot to aid in ventilation and cooling) and wait until the closure and carboy are completely cooled to room temperature before securing closure or closing spigot—**closure must be completely removed prior to autoclaving with no threads engaged**
- If autoclaving with liquid inside the carboy with spigots, remove spigot and replace with thread spigot closure (provided) in place over the threaded drain outlet—**do not autoclave with spigot in place as it will leak during autoclaving** (spigot can be autoclaved separately in an autoclave bag)
- If more than 1L of liquid is in the carboy, it may take several hours or more to reach the target sterilization temperature so it is usually more practical to autoclave the carboy and sterilize liquid by in-line filtration directly into the container
- Disposable carboy vent filter (Cat. No. 223-0030) is available featuring a hydrophobic, Teflon PTFE membrane, which permits sterile venting on slow exhaust/liquid autoclave cycles for the Nalgene PP and PC carboys up to 50L. This vent filter can be used up to 5 times

## Spigot Installation Instructions for Thermo Scientific Nalgene Carboys with Spigots

- Screw the lock nut of the spigot onto the threaded drain port of the carboy
- Lay the carboy on its back on a firm surface. Push the spigot down on the threaded drain port until you feel it snap into place and the lock nut again turns freely
- Adjust the spout so it points in the correct dispensing direction and turn the lock nut until it is tight again

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© 2013 Thermo Fisher Scientific Inc. All rights reserved. All (other) trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.