This product is molded of Teflon FEP (fluorinated ethylene propylene), a fluoropolymer which is inert to virtually all chemicals. The linerless, non-contaminating screw closure is molded of Tefzel ETFE (ethylenetetrafluoroethylene copolymer) and is designed to provide leakproof service. ETFE has similar chemical resistance to FEP and has even greater impact resistance and mechanical strength than FEP. Even under severe conditions this product is virtually unbreakable. Both FEP and ETFE are non-toxic and safe for contact with foods and beverages. They meet the requirements of the Food Additives Amendment of the U.S. Federal Food, Drug and Cosmetic Act.

Cleaning

Wash the labware before using for the first time. To prevent scratching of the plastic, observe the following precautions:

- Use a warm, mild detergent solution
- DO NOT use brushes, abrasive cleansers or paper towels
- If using a labware washing machine:
 - DO NOT use a machine that is equipped with brushes.
 - In a machine with high-pressure water spray place the labware in a basket and protect it with a screen or cover similar to that used in test tube baskets (such as Nalgene® autoclave baskets, Cat. No. 6917-0127, -0150 or -0230.). Otherwise, the water pressure may cause the labware to tumble and get scratched.

- If placed on spindles, always weight and cover the labware. If the spindles are made of uncoated metal, a section of PVC tubing placed over them will cushion and protect the labware.
- Air dry if possible.

Autoclaving

Clean the labware thoroughly before autoclaving. You can autoclave or sterilize FEP labware by any other thermal or chemical method.

NOTE: Before autoclaving, just set cap or closure on top of the container without engaging the threads.

Temperature Limits

This FEP container withstands temperatures from -270°C to +205°C. The ETFE closure withstands temperatures from -105°C to +150°C for extended periods. As with any plastic labware, **NEVER put in a flame or on a hot plate.**

For more information, see the current Nalgene Labware Catalog, or contact Technical Support.

© 2010 Thermo Fisher Scientific Inc. All rights reserved. Teflon and Tefzel are registered trademarks of DuPont. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

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