

Pierce<sup>®</sup> Glutathione Coated 96-Well Plates

15140 15240 15340

0685.2

Number	Description
15140	Pierce Glutathione Coated Plates (8-well strips), 5 plates
15240	Pierce Glutathione Coated Plates (white, 96-well), 5 plates
15340	Pierce Glutathione Coated Plates (black, 96-well), 5 plates

Capacity: ~10ng purified GST per well  
Blocking Buffer: Plates are supplied blocked with SuperBlock<sup>®</sup> Blocking Buffer  
Coating Volume: 100µL  
Blocking Volume: 200µL

**Storage:** Upon receipt store plates at 4°C in unopened pouches. Once opened, place unused plates in a resealable bag with desiccant and store at 4°C. Plates are shipped at ambient temperature.

## Introduction

Before purification of GST-fusion proteins using glutathione-agarose beads, it is often useful to screen for the presence of GST proteins in cell lysates. The Thermo Scientific Pierce Glutathione Coated Plates have glutathione immobilized on plates through its central sulfhydryl. These plates are useful for analyzing cell lysates to determine the presence and concentration of GST-fusion protein and for screening sera for antibodies to the GST-fusion protein. These plates are pre-blocked, making initial purification of the cell lysates unnecessary.

## General Procedure for GST-Fusion Protein Detection

### A. Materials Required

- Cell lysate containing GST-tagged protein
- Dilution buffer: Phosphate-buffered saline (PBS; Product No. 28372) or Tris-buffered saline (TBS; Product No. 28376)
- Wash Buffer: PBS or TBS with 0.05% Tween<sup>®</sup>-20 Detergent (Product No. 28320)
- Anti-GST antibody, 1mg/mL (Product No. 30001) or antibody against the protein of interest
- Enzyme- or fluorophore-conjugated secondary antibody to the species in which the anti-GST was created
- Appropriate enzyme substrate: example substrates are the Thermo Scientific TMB Substrate Kit (Product No. 34021) for horseradish peroxidase and the Thermo Scientific Phosphatase Substrate Kit (Product No. 37620) for alkaline phosphatase

### B. Method

1. Rinse each well three times with 200µL of wash buffer.
2. Prepare serial dilutions of the cell lysate in wash buffer.

**Note:** To determine the approximate concentration of GST-tagged protein, generate a standard curve using purified GST.

3. Apply 100µL of each lysate serial dilution to duplicate wells. Apply 100µL/well of wash buffer to separate duplicate wells to measure background. Cover the plate and incubate for 1 hour at room temperature.
4. Rinse the plate three times with 200µL/well of the wash buffer.

5. Prepare an appropriate dilution of the anti-GST antibody in dilution buffer and apply 100µL to each well. Cover the plate and incubate for 1 hour at room temperature.
6. Wash the plate three times with 200µL/well of wash buffer.
7. Prepare an appropriate dilution of the labeled secondary antibody and add 100µL to each well. Cover the plate and incubate for 1 hour at room temperature.
8. Wash the plate three times with 200µL/well of wash buffer.
9. Follow manufacturers' instructions for the specific detection system being used.

## Related Thermo Scientific Products

28372	BupH™ Phosphate Buffered Saline Packs, 40 packs
28376	BupH Tris Buffered Saline Packs, 40 packs
37070	SuperSignal® ELISA Pico Chemiluminescent Substrate, 100mL
15169	QuantaBlu™ Fluorogenic Peroxidase Substrate Kit
34028	1-Step Ultra TMB-ELISA, 250mL
37621	1-Step PNPP, 100mL
15075	Reagent Reservoirs, 200/pkg
15082	Microtube Racked System, 960 tubes
15036	Sealing Tape for 96-Well Plates, 100/pkg

## General References

- Smith, D.B. and Johnson, K.S. (1988). Single-step purification of polypeptides expressed in *Escherichia coli* as fusions with glutathione S-transferase. *Gene* 7:31-40.
- Frangioni, J.V. and Neel, B.G. (1993). Solubilization and purification of enzymatically active glutathione S-transferase (pGEX) fusion proteins. *Anal Biochem* 210:179-87.
- Simons, P.C. and VanderJagt, D.L. (1977). Purification of glutathione S-transferases for human liver by glutathione-affinity chromatography. *Anal Biochem* 82:334-41.

Tween is a trademark of Croda International Plc.

This product ("Product") is warranted to operate or perform substantially in conformance with published Product specifications in effect at the time of sale, as set forth in the Product documentation, specifications and/or accompanying package inserts ("Documentation") and to be free from defects in material and workmanship. Unless otherwise expressly authorized in writing, Products are supplied for research use only. No claim of suitability for use in applications regulated by FDA is made. The warranty provided herein is valid only when used by properly trained individuals. Unless otherwise stated in the Documentation, this warranty is limited to one year from date of shipment when the Product is subjected to normal, proper and intended usage. This warranty does not extend to anyone other than the original purchaser of the Product ("Buyer").

**No other warranties, express or implied, are granted, including without limitation, implied warranties of merchantability, fitness for any particular purpose, or non infringement. Buyer's exclusive remedy for non-conforming Products during the warranty period is limited to replacement of or refund for the non-conforming Product(s).**

There is no obligation to replace Products as the result of (i) accident, disaster or event of force majeure, (ii) misuse, fault or negligence of or by Buyer, (iii) use of the Products in a manner for which they were not designed, or (iv) improper storage and handling of the Products.

Current product instructions are available at [www.thermoscientific.com/pierce](http://www.thermoscientific.com/pierce). For a faxed copy, call 800-874-3723 or contact your local distributor.

© 2012 Thermo Fisher Scientific Inc. All rights reserved. Unless otherwise indicated, all trademarks are property of Thermo Fisher Scientific Inc. and its subsidiaries. Printed in the USA.