


Owner's Manual



Model D3-14



Model D2

The Model D3-14 Centipede™ & Model D2 Spider™ Horizontal Electrophoresis System

Manual No.:
D2_D3-0402
Rev. Date: 09/2002



Introduction

The Owl Model D3-14 Centipede™ & Model D2 Spider™ Horizontal Agarose Gel Electrophoresis Systems are designed to provide flat, even banding patterns and consistent results with hassle-free gel casting. No tape, grease, agarose seals or other accessories are required. A stand-alone casting platform is included for casting 1 or 2 (D2 only) gels simultaneously. Custom combs are available upon request.

Safety Information

Important Safety Information!
Please read carefully before operating!



- This manual contains important operating and safety information. In order to benefit from the use of this apparatus, you must carefully read and understand the contents of this manual prior to use of this apparatus.
- To avoid the risk of personal shock, always disconnect the gel box from the power supply. Further, the power supply must be equipped with a shut down on disconnect circuit.
- Statement of Proper Use: Use this product only for its intended purpose as described in this manual. Do not use this product if the power leads are damaged or if any of its surfaces are cracked.
- Running conditions for this unit should not exceed the name plate readings found on the lower buffer chamber.
- Do not move the unit unless the power source to the unit has been disconnected.



Table of Contents

Model D3-14 Centipede™ & D2 Spider™ Horizontal Electrophoresis System

Safety Information	i
Introduction	i
Unpack & Check Your Order.....	2
Model D2 Spider™	2
Model D3-14 Centipede™	3
Setting Up	4
Priming the Unit	4
Casting a Gel	4
Using the System	6
Running a Gel	6
Finishing Up	8
Care & Cleaning	9
Reagents Information	10
Selection of Reagents for Gel Electrophoresis	10
Troubleshooting	13
Specifications	15
Recommended Running Conditions	15
Migration Distance	15
Optional Equipment	16
Warranty Information	back cover

Unpack and Check Your Order


Before starting, unpack the unit and inventory your order. If any parts are missing, refer to the warranty section of this manual and contact Owl within 48 hours. Reference the order or catalog number on your invoice and check the corresponding part lists:

Model D2 Spider™ Wide Electrophoresis System

Gel Size: 14.4cmW x 10.2cmL
Footprint: 17cmW x 17cmL x 10cmH
Running Buffer Volume: 600ml

Complete System Includes:

- Buffer Chamber
- SuperSafe™ Lid with Attached Power Supply Leads
- EasyCast™ Gasketed U.V. Transmissible Gel Tray
- External Multiple Gel Caster
- 4 Combs: (2) 30 Well and (2) 40 Well, 1.5mm Thick



Item Description	Catalog No.
Complete System	D2
Complete System with Buffer Exchange Ports	D2-BP
Accessories	Catalog No.
Power Supply Leads	PSL-5
Gasketed EasyCast™ UVT Gel Tray	D2-UVT
Replacement Gaskets (1 pair)	D2-GK
External Multiple Caster (trays not included)	D2-CST
Leveling Platform 36cmW x 46cmL	B-LP
Bubble Level	BBL-1


Unpack and Check Your Order

Model D3-14 Centipede™ Wide Electrophoresis System

Gel Size: 23cmW x 14cmL
Footprint: 26cmW x 29cmL x 8cmH
Running Buffer Volume: 800ml

Complete System Includes:

- Buffer Chamber
- SuperSafe™ Lid with Attached Power Supply Leads
- EasyCast™ Gasketed U.V. Transmissible (UVT) Gel Tray
- External Gel Caster
- 4 Combs: 50 Well 1.5mm Thick



Item Description	Catalog No.
Complete System	D3-14
Accessories	
Power Supply Leads	PSL-5
EasyCast™ Gasketed UVT Gel Tray	D3-UVT-14
Replacement Gaskets (1 pair)	D3-GK
External Gel Caster (trays not included)	D3-CST-14
Leveling Platform 36cmW x 46cmL	B-LP
Bubble Level	BBL-1

Setting Up

Item	Part No.	Description	Qty
Buffer chamber		Safety lid with attached power cords	
Gasketed UVT gel tray			

Casting chamber
4 combs, 1.5mm, 2 each 30 and 40 wells

Casting a Gel

STEP 1

Remove the SuperSafe lid from the gel box by holding the front of the buffer chamber with one hand and pulling the lid off; sliding off evenly by holding the center of the back of the lid. The cover is attached to the back of the unit at the junction of the lids attached power cords to the banana plugs located on the unit.



D2

D3-14

STEP 2

For shipping and convenient storage, the gasketed gel tray is packaged inside the casting chamber upon arrival.

To remove the gel tray, lift the tray out of the



caster by placing both hands firmly on the shorter tray ends and pull up slowly

from the caster at an angle. Tray needs to fit snug for leakproof gel casting, so it may be tight. "Walking" the tray upward at an angle may be helpful.

STEP 3

To cast gels, place the gel tray into the external gel caster making sure the tray is pushed all the way down and level in the external gel caster. For best results, be sure to cast on a level surface. Leveling platforms are available if required, Catalog #B-LP.

STEP 4

Setting Up





The need to cast a level gel is very important for consistent reproducible results. Level the caster by placing the enclosed bubble level (Model D3-14 only, order B-LB for Model D2) in the center of the gel tray. Using the thumbscrews on each side of the caster, slowly turning one thumbscrew (front only) at a time and lining up the bubble in the level with the center circle. Check various areas on the gel tray by moving the bubble level to each end of the gel tray to ensure you have leveled correctly.

NOTE: It is wise to always run a sample lane of a known “standard ladder” to determine concentration and size of separated fragments after the gel run, and to aid in photo documentation and analysis.

STEP 5

Preparing the gel - Using electrophoresis grade agarose and compatible electrophoresis buffer the gel may be prepared in various ways. The percentage of agarose and the electrophoretic buffer used is determined by the size of the samples to be separated and further recovery of the samples (see page 9 and 10). The agarose and buffer are mixed and heated using a heat source, in a microwave oven, or in an autoclave until the agarose is completely dissolved. The

prepared gel then must be cooled to below 60° before casting to avoid warping the UVT gel tray due to excessive heat. If numerous gels are to be run in one day, a large volume of gel may be prepared and placed in a covered bottle stored between 40-60° in a water bath. This provides a ready gel supply in a warm liquid form that will solidify quickly when gels are cast. For further tips on sample preparation and visualization, see page 9 and 10.



STEP 6

Pour or pipet the correct amount (see page 9 and 10) of warm agarose (< 60°) onto the UVT gel tray. Immediately after pouring, insert the desired comb or

Setting Up

combs into the comb slots to form the sample wells. If only a small portion of gel is required for proper sample separation, up to 4 combs may be used to run 2, 3 or 4 (see page 14) sets of equal distance samples simultaneously expanding the number of samples per gel that may be run. Note: A higher percentage of agarose (>0.5%) may give the best results when using 4 of the 50 tooth (D3-14 only, D2 goes to 40) combs to avoid damage to the wells when combs are removed.

Wall Comb: To conserve agarose, a wall comb may also be used to divide and use a shorter length of the gel tray. If a wall comb is used, pipet a bead of agarose along the bottom and side edges of the wall comb once it has been placed in the tray to seal the

Using the System

Running A Gel

STEP 1

Once the gel is completely solidified, carefully lift the tray (2) out of the external gel caster as described above and place into the buffer chamber. The running position of the tray exposes the open ends of the agarose to the buffer.

STEP 2

Pour enough compatible running buffer into the unit to fill both ends of buffer chamber and completely cover and submerge the gel. Correct buffer level is clearly marked on the side wall as "FILL LINE". See page 14, Recommended Running Conditions, for approximate buffer volumes needed for your unit. Too little buffer may cause the gel to dry out during the run, while excess buffer may decrease DNA mobility and cause band distortion.

Using the System

STEP 3

Carefully remove the comb (or combs) by tapping lightly to loosen, and slowly lifting straight up out of the gel tray. To avoid damage to the sample wells, always make sure to allow the gel to solidify completely before moving the buffer chamber, gel tray, or removing the combs. After placing the gel tray into the unit in the running position, submerge the gel in 3-5mm of running buffer. Lightly tap each comb gently back and forth to loosen, then slowly pull the comb straight up out of the gel tray. This will break any suction that may exist between the gel and comb. When using all four combs (D3-14 model), a higher percentage of agarose (>0.5%) may be wise to avoid damage to the sample wells. Low percentage gel and the small sample wells may cause the sides of the wells to collapse when the comb is removed. A higher percentage of agarose forms a tighter gel matrix. Casting a slightly thicker gel may also remedy this problem.

NOTE: Combs may also be removed prior to pouring buffer in buffer chamber.

STEP 4

Load prepared samples into the wells (see page 14). Samples should be mixed with a sample loading buffer; giving weight to the samples so that they drop evenly into the wells and contain tracking dyes to monitor the gel run. Refer to page 16. The combs supplied with the D3-14 unit are designed in a micro well format. This format allows speedy sample loading using a multi-channel pipet. The 25 tooth comb is in the 1X micro well format and matches each tip of the multichannel pipet, while the 50 tooth comb is in the 2X format; loading every other lane. The D2 has a 1X (15 well comb) and a 2X (30 well comb).

STEP 5

Carefully slide the SuperSafe lid with attached power supply leads onto the unit. This will connect the power supply leads to the banana plug electrodes and complete the circuit. Plug other end of the power cords into appropriate power supply.

STEP 6

Turn on the power supply (See page 14, Recommended Running Conditions). Carefully monitor the gel run to avoid samples running into the path of another set of samples.

Finishing Up

1. When the gel run is complete and tracking dye has migrated as far through the gel as desired or to the end of the gel, turn off the power supply and slide off the super safe lid to disconnect from the power source. Carefully remove the tray containing the gel (wear gloves if ethidium bromide is present). The UV transparent gel tray makes visualization and photography with a UV light source easy without the need to remove the gel from the tray. The gel tray may be placed back into the casting chamber for convenient transport to the darkroom to avoid damage to the gel.
2. The gel box should be rinsed under warm running water after each use, including the UVT gel tray. Rinsing the gel tray will avoid any salt build up in the gasket material from the electrophoretic running buffer; extending the gasket life and ensure leak-free gel casting.

Care & Cleaning

A Few Tips About Caring for Your System

WARNING!

DO NOT USE ETHANOL OR OTHER ORGANIC SOLVENTS TO CLEAN OWL PRODUCTS. Organic solvents cause acrylic to "craze" or crack. Clean all Owl acrylic systems with warm water and a mild detergent.

Do not autoclave, bake, or microwave your unit. Temperatures over 50°C can do damage to the acrylic.

The unit may be rinsed with warm water, or cleaned with warm water and a mild detergent to get rid of any debris.

NOTE:

If an RNase free electrophoresis system is desired, there are various methods to rid the system of RNA contamination. For fast and easy decontamination, use RNase Away®*. Spray, wipe or soak labware with RNase Away® then wipe or rinse the surface clean; it instantly eliminates RNase. RNase Away® eliminates the old methods that include treatment with 0.1% Diethyl Pyrocarbonate (DEPC) treated water and soaking in dilute bleach. DEPC is suspected to be a carcinogen and should be handled with care. This electrophoresis system should never be autoclaved, baked, or placed in a microwave.

To order RNase Away®, contact Molecular BioProducts 800-995-2787 (U.S. and Canada) or 858-453-7551:

Part Number

7000	250ml bottle
7002	475ml spray bottle
7003	1 liter bottle
7005	4 liter bottle

*Rnase AWAY® is a registered trademark of Molecular BioProducts

Reagents Information

SELECTION OF REAGENTS FOR GEL ELECTROPHORESIS

1. Agarose

There are various types of agarose commercially available that may be used. Besides standard ultra pure electrophoresis grade agarose, there are also numerous low melting point products for easy sample recovery, as well as specialty products formulated for specific uses (i.e. to separate and/or recover very small or very large fragments).

Table 2: Mobility range of DNA in different percentage agarose gels

Agarose % (w/v)	Approximate range of separated DNA fragments (kb)
0.3	60 to 5
0.5	30 to 1
0.7	12 to 0.8
1.0	10 to 0.5
1.2	7 to 0.3
1.5	4 to 0.2
2.0	3 to 0.1
3.0	<0.1

Table 3: Amount of Agarose to prepare

Gel volume is determined by the following formula and may be adjusted according to need or preference: $gel\ width(cm) \times gel\ length\ (cm) \times gel\ thickness\ (cm) = ml\ of\ agarose$

Model #	Gel size(cm)	0.25cm	0.5cm	0.75cm	1.0cm
D3-14	23x14	80.5ml	161ml	241.5ml	322ml
D2	14.4x10.2	36.72	73.44	110.16	146.88

NOTE: an increased agarose % provides better separation of small fragments and bands very close together that tend to be more difficult to separate. A specialty agarose product formulated to increase resolution of low molecular mass samples may also be used, or an agarose additive may be added to standard or low melting point agarose.

Example: A good mid range gel percentage would be 0.7%, or 0.7g agarose in 100ml electrophoresis buffer (TBE or TAE), following heating and dissolving the agarose, 10ul. of ethidium bromide stock solution (5mg/ml) is added. The gel would be run with compatible electrophoretic running buffer (1X TBE or 1X TAE) that also contained ethidium bromide. One liter of the running buffer would contain 100ul. of this 5mg/ml. ethidium bromide stock solution.

Reagents Information

2. Ethidium Bromide

For photodocumentation of samples, the gel may be stained during or following the run with a variety of stains. The most common stain for DNA is ethidium bromide. Ethidium bromide may be added directly to the gel and running buffer to visualize and photograph the separated fragments following the gel run without the need for an additional staining step. The ethidium bromide is added to both the gel (after heating) and the electrophoresis buffer at a concentration of 0.5ug/ml. Conversely, the gel may be stained in a concentrated ethidium bromide solution after the gel run and rinsed for visualization.

NOTE: Ethidium bromide is a potential carcinogen. Care in handling the powder and stock solution must be taken. Always wear gloves when handling the powder, solutions and all gels that contain ethidium bromide.

**Table 5: Preparation & Properties of
TAE and TBE Electrophoresis Buffer Systems**

These buffers are used because they both have a basic pH which gives the phosphate group of the DNA a net negative charge allowing migration of the DNA toward the positive anode in the electrophoresis chamber.

TAE - Tris Acetate with EDTA (40mM Tris Base, 40mM Acetic Acid, 1mM EDTA)

<u>50X stock solution, pH 8.5</u>	<u>1X working solution:</u>
242g Tris Base	40mM Tris Acetate
57.1ml Glacial Acetic Acid	1mM EDTA
18.61g Na EDTA " 2H ₂ O (MW 372.24)	
Distilled H ₂ O to 1 Liter Final Volume	

TBE - Tris Borate with EDTA (89mM Tris Base, 89mM Boric Acid, 2mM EDTA)

<u>10X stock solution:</u>	<u>1X working solution:</u>
108g Tris Base	89mM Tris Base
55g Boric Acid	89mM Boric Acid
7.44g Na ₂ EDTA " 2H ₂ O (MW 372.24)	2mM EDTA
(or 40ml 0.5M EDTA, pH 8.0)	
Distilled H ₂ O to 1 Liter Final Volume	

Reagents Information

Choose the buffer best suited to the experiment. Each buffer has different properties providing the necessary ions for electrophoretic migration.

Buffer:

Suggested Use:

TAE Buffer

- Use when DNA is to be recovered
- For electrophoresis of large (>20kb) DNA
- Applications requiring high resolution
- Has low ionic strength and low buffering capacity - recirculation may be necessary for long runs (>4hrs.)

TBE Buffer

- General Purpose Buffer
- Can be re-used
- For electrophoresis of small (<1kb) DNA
- Better resolution of small (<1kb) DNA
- Decreased DNA mobility
- High ionic strength and high buffering capacity - recirculation may not be required for extended run times
- Reacts with the agarose making smaller pores and a tighter matrix. This reduces broadening of the DNA bands for sharper resolution.

3. Sample Buffer

Samples are prepared and mixed with sample buffer before being applied to the prepared gel. Sample buffers contain similar components to the running buffer, dyes for visibility, and glycerol to provide weight to the samples. This increased sample density ensures samples load evenly into the wells and do not float out during loading. Dyes also migrate toward the anode end of the electrophoresis chamber at predictable rates allowing the gel run to be monitored.

4. DNA Markers

Markers are run on each gel to monitor sample separation and to provide an accurate size estimation of the samples. By running a known marker of a specific concentration, the amount of the DNA can be estimated. These size markers are a suitable restriction digest of commonly available DNA.

Troubleshooting

PROBLEM

SOLUTION

Agarose leaks into chamber when pouring gel

- Check to see if the gasket is firmly seated in the grooves on the ends of the UVT gel tray. Reseat gasket if necessary by removing and rinsing under warm running water, then reseat evenly in the tray groove.

Bands seem to be running at an angle.

- Check to be sure the casting is being done on a level surface. A leveling platform may be required. Make sure the UVT gel tray is pressed all the way down and rests level on the external gel casting platform (the bubble in the bubble level should rest in the center circle). Adjust the leveling screws to make the external gel caster level.

Samples seem to be running unevenly in certain areas.

- Check to be sure the platinum electrode wire is intact and running evenly across the base of the chamber and up the side to the junction of the banana plug. If there appears to be a break in the electrode connection contact Owl immediately. This problem may also be caused by regular casting with very hot agarose gel (>60°F) which may damage the UVT gel tray over time. Always cool the melted agarose to below 60°F before casting to avoid warping the UVT gel tray. Warping the gel tray will cause all subsequent gels to be cast unevenly.

Samples do not band sharply and appear diffuse in the gel.

- Gels should be no more than 5mm thick and allowed to solidify completely before running. For standard agarose this would be about 30 minutes, if low melting point agarose is used it may be necessary to completely solidify gels at a cooler temperature in the refrigerator or cold room. Gels should be submerged in 3-5mm of buffer to avoid gel dry out, but excess buffer >5mm can cause decreased DNA mobility and band distortion.

Samples are not moving as expected through the gel, remaining in the wells, running "backwards" or diffusing into the gel.

- Check to be sure that a complete power circuit is achieved between the unit and the power supply. Platinum wire and banana plugs should be intact. To test, simply fill the buffer chamber with running buffer and attach to the power supply leads without a gel or UVT gel tray in the unit. The platinum wires on both sides of the buffer chamber should produce small bubbles as the current passes through. If a complete circuit does not exist there will be little to no bubbles. Contact Owl's Customer Service Department to schedule a repair. Samples that appear to run backwards through the gel is caused by the UVT tray being placed in the buffer chamber in the reverse direction. The tray should be placed in the buffer chamber with the comb at the edge of the tray closest to the cathode side of the buffer chamber.

Troubleshooting

PROBLEM

SOLUTION

When the comb is removed from the gel the sample well is ripped and damaged.

- Always make sure to allow the gel to solidify completely before moving the UVT tray, buffer chamber or removing the comb. To avoid damage to the sample wells, gently rock the comb back and forth lightly to loosen, then slowly pull the comb straight up out of the UVT gel tray. This rocking helps to avoid suction as the comb is removed.

The gel seems to run slower under the usual running conditions.

- The volume of running buffer used to submerge the gel should only be between 3-5mm over the gel surface. The gel should be completely submerged to avoid the gel from drying out, which can smear the bands and possibly melt the gel due to overheating. If excessive running buffer is added the mobility of the DNA decreases and band distortion may result. Excess buffer causes heat to build up and buffer condensation inside the unit may result.

Additional Sources for Reference

Maniatis T., E. F. Fritsch and J. Sambrook. *Molecular Cloning: A Laboratory Manual*, Cold Spring Harbor Laboratory, Cold Spring Harbor, NY.

Short Protocols in Molecular Biology, - A Compendium of Methods from Current Protocols in Molecular Biology, Edited by Fredrick M. Ausubel, et. al.

Adams, D., and R. Ogden, *Electrophoresis in Agarose and Acrylamide Gels, Methods in Enzymology*, Vol. 152 (1987) Academic Press, Inc.

Fotador, U.. *Simultaneous Use of Standard and Low-Melting Agarose for the Separation and Isolation of DNA by Electrophoresis*, BioTechniques, Vol. 10, No. 2, (1991)

Boots, S. *Gel Electrophoresis of DNA* ; Analytical Chemistry, Vol. 61, No. 8, April 15, 1989

Specifications





RECOMMENDED RUNNING CONDITIONS

Model	D3-14	D2
Gel Size (W x L in cm)	23x14	14.4 x 10.2
Buffer Capacity (L)	800ml	600ml
Voltage Requirements (V)	20-150	20-150
Time Requirements (min)	30-60	30-60

Migration Distance


The charts below give the migration distance for each comb slot on the gel tray with the run lengths.

D3-14 Migration Distance




Four combs for four 3cm run lengths

D2 Migration Distance



One comb for a 9.2cm run length




Two combs for two 4.4cm run

Optional Equipment

Buffer Exchange Port Option, D2-BP (Model D2 only)

The buffer exchange port option is used to recirculate the buffer during extended gel runs. Recirculation is used to prevent buffer depletion of certain low ionic running buffers, for extended runs multiple sample sets, or for RNA gels. If your unit has the buffer exchange port option it will be fitted with two white buffer port terminals and will contain two separate port inserts packaged in a small plastic bag located inside the unit upon arrival.



How these work...

The inserts are pushed into the attached ports on the side wall of the unit with the black O-ring side facing in. The insert will “snap” into place in the port in the “open” position and is ready to circulate buffer. Appropriate tubing is then connected to the small outer ringed ends of the ports for circulation using a separate recirculator or peristaltic pump. To close the port, which also releases the insert, you simply press the flat metal button and the insert detaches. The port is now in the “closed” position.

NOTE: Buffer may also be passed through a heat exchanger.

Comb Options

The following charts list the available comb options for the D2 and D3-14.

Comb Options – Model D2					Recommended Loading Volumes ¹			
Catalog Number	Comb Type	Number of Teeth	Thickness of Teeth	Width of Teeth	0.25cm ²	0.5cm ²	0.75cm ²	1.0cm ²
D1-10C	Standard	10	1.0 mm	12.2mm	9 ul	32 ul	55 ul	78 ul
D1-10D	Standard	10	1.5	12.2	14	48	82	117
D1-20C	Standard	20	1.0	5.2	4	14	23	33
D1-20D	Standard	20	1.5	5.2	6	20	35	50
D1-30C	Micro Well	30 (2x) ³	1.0	7.2	5	19	32	46
D1-30D	Micro Well	30 (2x) ³	1.5	7.2	8	28	49	69
D1-40C	Standard	40	1.0	1.7	1	4	8	11
D1-40D	Standard	40	1.5	1.7	2	7	11	16
D1-MTC	Micro Well	15 (1x) ³	1.0	7.2	5	19	32	46
D1-MTD	Micro Well	15(1X) ³	1.5	7.2	8	28	49	69
XCM	Custom		1.0, 1.5 2.0, 3.0					

¹ Loading Volume is calculated as 75% of total well volume (TxWxHx0.75)

² Gel Thickness

³ 8 & 12 Channel Pipette Format

Optional Equipment

Comb Options – Model D3-14					Recommended Loading Volumes ¹			
Catalog Number	Comb Type	Number of Teeth	Thickness of Teeth	Width of Teeth	0.25cm ²	0.5cm ²	0.75cm ²	1.0cm ²
					D3-MTC	Micro Well	25 (1X) ³	1.0mm
D3-MTD	Micro Well	25 (1X) ³	1.5	7.5	8	30	51	72
D3-MT2C	Micro Well	50 (2X) ³	1.0	3	2	8	14	19
D3-MT2D	Micro Well	50 (2X) ³	1.5	3	3	12	20	29
D3-WALL	Wall	1		230				
XCM	Custom		1.0, 1.5 2.0, 3.0					

¹ Loading Volume is calculated as 75% of total well volume (TxWxHx0.75)

² Gel Thickness

³ 8 & 12 Channel Pipette Format

How to determine well sample volume:


Hg = height of gel used Hs = height of well used for sample volume Hw = well height

There are two volumes to consider when determining the sample volume for a horizontal gel.

1) Gel volume, which is Width x Length x Gel Height and uses centimeters and 2) Sample volume which is Tooth Width x Comb Thickness x Apparent Well Height, and uses millimeters.

Gel Height is generally set to a height between 0.5 cm and 0.75 cm. Therefore, once you choose the height, the volume is the gel dimensions (given in the catalog for each gel box, I.D.) times this height.

Once the gel height (Hg) is chosen, the well volume and then the sample volume can be calculated.



The well height (Hw) is 1.5 mm less than the gel height: $Hw = \text{Gel Height} - 1.5 \text{ mm}$

Using the well height, the volume of the well is calculated: $Vw = (\text{Well Height}) \times (\text{Tooth width} \times \text{comb thickness})$

The loading volume is a 0.75 safety factor applied to the well volume: $Vs = (Vw) \times (.75)$

For Owl Combs, there are only two thicknesses, 1.0mm and 1.5mm. This is the depth. The width of the well is determined by the number of teeth. For a given gel box, as the number of teeth increase, the volume of each tooth decreases.

Warranty Information

The Owl Separation Systems Warranty

A three-year quality and material warranty covers all products manufactured by Owl Separation Systems. Owl will repair or replace any equipment found to be defective at no cost. This warranty does not cover equipment damage due to misuse or abuse. After the warranty expires, Owl will repair products at a reasonable cost. All shipping claims must be made within 48 hours from date received.

To activate your warranty, complete and return the enclosed postage paid warranty card. Please note that the card must be completely filled out in order to process your warranty.

Returning Equipment

Be environmentally friendly – and speed up your return – by saving all packing materials cartons and documents until you have thoroughly inspected your shipment. Should you find that your order is incorrect or damaged, verify the problem with the shipper, save all packing material, and call Owl for return instructions within 48 hours. All returns, exchanges, and credits must be pre-approved by Owl.



55 Heritage Avenue
Portsmouth, NH 03801

T. (603) 559-9297
(800) 242-5560
F. (603) 559-9258

Website: www.owlsci.com
E-mail: sales@owlsci.com

Thank You!

We at Owl Separation Systems thank you for your order and appreciate your business. Please contact us regarding our complete line of electrophoresis equipment and reagents for DNA, RNA and protein separations. While innovation and quality are our foremost objectives, we pride ourselves on exceptional customer response and service.