


CATALOG NUMBER	CLOSURE SIZE	MATERIALS				
		CLOSURE	GASKET	TUBING	PORT CAP	VENT CAP
2162-0531	53B (1/4" PORT)	PP (WHITE)	TPE	NALGENE 50 PLATINUM CURED SILICONE	TPE	TPE

THREAD FINISH (mm).....:	53	
THREAD (T.P.I.).....:	4	
NOMINAL WEIGHT (ASSEMBLED gm).....:	62	
	INCH	MM
A (NOMINAL HEIGHT).....:	2.7	69
D (DIA.).....:	2.59	65.5
BARBED PORT SIZE (INCH)....:	1/4	
VENT PORT SIZE (INCH).....:	1/4	

**NOTE:**

53B CLOSURES ARE SUPPLIED WITH A 4" SHORT TUBE AND A 18" LONG TUBE ATTACHED.

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



SPECS ON BACK

NOT TO SCALE

**CONFIDENTIAL**

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

**ThermoFisher**  
SCIENTIFIC

Thermo Fisher Scientific  
75 Panorama Creek Drive  
Rochester, NY 14625-2385, USA

MADE	PB	APP'D	GRP	PSS 00-2162-0531	REV.	0413
CK'D	KOW	4/14/13				

DESCRIPTION: 53B FILLING/VENTING CLOSURE