


CATALOG NO.	NOMINAL CAPACITY (GAL)	MATERIALS	
		JERRICAN	CLOSURE
002240	1.5	HDPE	WHITE PP

SIZE CODE.....:	-0015	
SIZE (NOM VOL L/GAL):	6/1.5	
NOM. BRIM CAP. (L)...:	7.6	
NECK FINISH (mm)....:	53B	
THREAD (T.P.I.).....:	4	
NOM. WEIGHT W/ CLOSURE(gms).....:	920	
	INCH	MM
A (NECK I.D.).....:	1.41 ±0.04	35.8 ±1.0
B (HT. W/ CLOSURE)...:	13.2 ±0.2	335 ±5
C (HT. W/O CLOSURE)..:	12.9 ±0.2	328 ±5
D (WIDTH).....:	8.38 ±0.14	212.9 ±3.6
E (DEPTH).....:	6.94 ±0.12	176.3 ±3.0
W.T.(NOM WALL THK)...:	0.18	4.6

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



SPECS ON BACK

NOT TO SCALE

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

**ThermoFisher**  
SCIENTIFIC

Thermo Fisher Scientific  
75 Panorama Creek Drive  
Rochester, NY 14625-2385, USA

MADE	TL	APP'D	GRP	PSS 00-2240-0015-S	REV.
CK'D	KOW		5/15/13		0413
DESCRIPTION:					
JERRICANS					