


CATALOG NO.	NOMINAL CAPACITY (L)	MATERIALS		
		CARBOY	CLOSURE	STRAP
002243-0013	13	HDPE	WHITE PP/NAT	HDPE

SIZE (NOM VOL L).....	13	
NOM BRIM CAP (L).....	14.2	
NECK FINISH (mm).....	53	
THREAD (T.P.I.).....	4	
NOM. WEIGHT W/ CLOSURE(gms).....	1250	
	INCHES	MM
A (NECK I.D. AT TOP)..:	1.5	38
B (OVERALL HEIGHT)..:	14.9	378
C (WIDTH).....	11.4	290
D (DEPTH).....	7.45	189
W.T.(NOM WALL THK)..:	0.1	2.5

2. This is not a production drawing, but is for informational purposes only, and is subject to change without notice.

1. All dimensions are nominal and in inches and millimeters unless otherwise noted.



SPECS ON BACK

NOT TO SCALE

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

**ThermoFisher**  
SCIENTIFIC

Thermo Fisher Scientific  
75 Panorama Creek Drive  
Rochester, NY 14625-2385, USA

MADE	TL	APP'D	GRP	PSS 00-2243-0013	REV. 0413
CK'D	KOW		5/15/13		

DESCRIPTION:  
13-LITER JERRICAN