


CATALOG NO.....:	6177-0250	6177-0375		
NOM. TUBE I.D.....:	1/4"	3/4"		
MATERIALS	BODY.....:	PP		
	SLEEVE.....:	DELIN (ACETAL)		
	SPRINGS & LATCHES...:	STAINLESS STEEL		
COLOR.....:	ALMOND			
	INCH	MM	INCH	MM
A (LENGTH).....:	3.28	83.3	3.09	78.5
B (BARB O.D.).....:	0.30	7.6	0.42	10.7



SPECS ON BACK

NOT TO SCALE

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

ThermoFisher
SCIENTIFIC

Thermo Fisher Scientific
75 Panorama Creek Drive
Rochester, NY 14625-2385, USA

MADE	MC	APP'D	PVC	PSS 00-6177-0XXX	REV.	1099
CK'D	PVC		10/13/99			

DESCRIPTION:

QUICK DISCONNECT BARBED VALVE