

CATALOG NUMBER	NOMINAL SIZE (mL)								MATERIAL
1510	30	50	100	150	250*	400*	600*	1000*	PFA



(*) FOR SIZES: -250, -400, -600, -1000 SEE PSS 00-1510-XXXX-2

SIZE CODE.....	-0030		-0050	
SIZE (NOM.VOL.mL)...	30		50	
NOM. BRIM CAP (mL)...	44		76	
NOM. WT. (gm).....	27		38	
	INCH	MM	INCH	MM
A (LIP DIA.).....	1.77 ±0.03	44.9 ±0.8	2.06 ±0.03	52.3 ±0.8
B (DIM.).....	1.88 ±0.03	47.7 ±0.8	2.24 ±0.03	56.9 ±0.8
C (I.D. TOP).....	1.43 ±0.02	36.3 ±0.5	1.71 ±0.03	43.4 ±0.8
D DIA.....	1.43 ±0.02	36.3 ±0.5	1.67 ±0.03	42.4 ±0.8
H (HEIGHT).....	1.95 ±0.03	49.5 ±0.8	2.32 ±0.03	58.9 ±0.8
W.T.(NOM. WALL THK):	0.07	1.8	0.07	1.8

SIZE CODE.....	-0100		-0150	
SIZE (NOM.VOL.mL)...	100		150	
NOM. BRIM CAP (mL)...	140		224	
NOM. WT. (gm).....	60		90	
	INCH	MM	INCH	MM
A (LIP DIA.).....	2.53 ±0.04	64.3 ±1.0	2.85 ±0.04	72.4 ±1.0
B (DIM.).....	2.72 ±0.04	69.1 ±1.0	3.13 ±0.05	79.5 ±1.3
C (I.D. TOP).....	2.10 ±0.03	53.3 ±0.8	2.46 ±0.04	62.5 ±1.0
D DIA.....	2.04 ±0.03	51.8 ±0.8	2.35 ±0.04	59.7 ±1.0
H (HEIGHT).....	2.80 ±0.04	71.1 ±1.0	3.24 ±0.05	82.3 ±1.3
W.T.(NOM. WALL THK):	0.07	1.8	0.09	2.3

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



SPECS ON BACK

NOT TO SCALE

NOT TO SCALE

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

ThermoFisher
SCIENTIFIC

Thermo Fisher Scientific
75 Panorama Creek Drive
Rochester, NY 14625-2385, USA

MADE	MG	APP'D	PVC	PSS 00-1510-XXXX-1	REV. 0695
CHK'D	PVC		6/15/95		

DESCRIPTION:

GRIFFIN LOW-FORM BEAKERS