

CATALOG NUMBER	NOMINAL SIZE (mL)								MATERIAL
1510	30 *	50 *	100 *	150 *	250	400	600	1000	PFA



(*) FOR SIZES: -30, -50, -100, -150 SEE PSS 00-1510-XXXX-1

CATALOG NO.....:	-0250		-0400	
SIZE (NOM.VOL.mL):	250		400	
NOM. BRIM CAP (mL):	340		500	
NOM. WT. (gm):	140		180	
	INCH	MM	INCH	MM
A (LIP DIA.):	3.28 ±0.05	83.3 ±1.3	3.77 ±0.06	95.8 ±1.5
B (DIM.):	3.55 ±0.05	90.2 ±1.3	4.05 ±0.06	102.9 ±1.5
C (I.D. TOP):	2.90 ±0.05	73.7 ±1.3	3.37 ±0.05	93.2 ±1.3
D DIA.....:	2.75 ±0.04	69.9 ±1.0	3.16 ±0.05	80.3 ±1.3
H (HEIGHT):	3.84 ±0.06	97.5 ±1.5	4.19 ±0.06	106.4 ±1.5
W.T.(NOM. WALL THK):	0.10	2.5	0.10	2.5

CATALOG NO.....:	-0600		-1000	
SIZE (NOM.VOL.mL):	600		1000	
NOM. BRIM CAP (mL):	830		1430	
NOM. WT. (gm):	265		342	
	INCH	MM	INCH	MM
A (LIP DIA.):	4.42 ±0.06	112.3 ±1.5	5.26 ±0.08	133.6 ±2.0
B (DIM.):	4.75 ±0.07	120.7 ±1.8	5.62 ±0.08	142.7 ±2.0
C (I.D. TOP):	3.95 ±0.06	100.3 ±1.5	4.61 ±0.07	117.1 ±1.8
D DIA.....:	3.62 ±0.06	91.9 ±1.5	4.30 ±0.06	109.2 ±1.5
H (HEIGHT):	5.02 ±0.08	127.5 ±2.0	5.90 ±0.08	259.1 ±2.0
W.T.(NOM. WALL THK):	0.11	2.8	0.10	2.5

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



SPECS ON BACK

NOT TO SCALE

NOT TO SCALE

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

ThermoFisher
SCIENTIFIC

Thermo Fisher Scientific
75 Panorama Creek Drive
Rochester, NY 14625-2385, USA

MADE	MG	APP'D	PVC	PSS 00-1510-XXXX-2	REV. 0695
CHK'D	PVC		6/15/95		

DESCRIPTION:

GRIFFIN LOW-FORM BEAKERS