


CATALOG NO.	NOMINAL CAPACITY, (oz)							MATERIALS	
								BOTTLE	CLOSURE
002019	1	2	4	8	16	32	64	PETG	WHITE HDPE

SIZE CODE.....	-0030	-0060	-0125	-0250				
SIZE (NOM VOL mL/oz)	30/1	60/2	125/4	250/8				
NOM. BRIM CAP. (mL):	45	77	175	320				
NECK FINISH (mm)...	20-415	24-415	38-430	38-430				
THREAD (T.P.I.).....	8	8	6	6				
NOM. WEIGHT W/CLOSURE (gms).....	14.5	22	52.5	62.8				
	INCH	MM	INCH	MM	INCH	MM	INCH	MM
A (NECK I.D.).....	0.54 ±0.02	13.7 ±0.5	0.70 ±0.02	17.8 ±0.5	1.100 ±0.025	27.94 ±0.64	1.100 ±0.025	27.94 ±0.64
B (HT. W/CLOSURE)...	2.50 ±0.05	63.5 ±1.3	3.25 ±0.06	82.6 ±1.5	4.26 ±0.08	108.2 ±2.0	5.70 ±0.08	144.8 ±2.0
C (HT. W/O CLOSURE):	2.39 ±0.05	60.7 ±1.3	3.15 ±0.06	80.0 ±1.5	4.08 ±0.08	103.6 ±2.0	5.50 ±0.08	139.7 ±2.0
D (SQUARE).....	1.50 ±0.05	38.1 ±1.3	1.58 ±0.05	40.1 ±1.3	2.06 ±0.06	52.3 ±1.5	2.32 ±0.06	58.9 ±1.5
R (CORNER RADIUS)...	0.21 ±0.03	5.3 ±0.8	0.21 ±0.03	5.3 ±0.8	0.22 ±0.04	5.6 ±1.0	0.31 ±0.04	7.9 ±1.0
W.T.(NOM. WALL THK):	0.04	1.0	0.04	1.0	0.06	1.5	0.06	1.5

SIZE CODE.....	-0500	-1000	-2000			
SIZE (NOM VOL mL/oz)	500/16	1000/32	2000/64			
NOM. BRIM CAP. (mL):	610	1225	2380			
NECK FINISH (mm)...	38-430	38-430	53B			
THREAD (T.P.I.).....	6	6	4			
NOM. WEIGHT W/CLOSURE (gms).....	116	187	350			
	INCH	MM	INCH	MM	INCH	MM
A (NECK I.D.).....	1.100 ±0.025	27.94 ±0.64	1.100 ±0.025	27.94 ±0.64	1.53 ±0.03	38.9 ±0.8
B (HT. W/CLOSURE)...	7.00 ±0.09	177.8 ±2.3	8.63 ±0.10	219.2 ±2.5	10.65 ±0.14	270.5 ±3.6
C (HT. W/O CLOSURE):	6.82 ±0.09	173.2 ±2.3	8.48 ±0.10	215.4 ±2.5	10.42 ±0.14	264.7 ±3.6
D (SQUARE).....	2.88 ±0.06	73.2 ±1.5	3.63 ±0.07	92.2 ±1.8	4.55 ±0.09	115.6 ±2.3
R (CORNER RADIUS)...	0.38 ±0.04	9.7 ±1.0	0.50 ±0.04	12.7 ±1.0	0.75 ±0.05	19.1 ±1.3
W.T.(NOM. WALL THK):	0.065	1.7	0.05	1.3	0.09	2.3

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



SPECS ON BACK

NOT TO SCALE

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

ThermoFisher
SCIENTIFIC

Thermo Fisher Scientific
75 Panorama Creek Drive
Rochester, NY 14625-2385, USA

MADE PB	APP'D PWC	PSS 00-2019-XXXX-1	REV. 0514
CK'D KOW	5/20/14		
DESCRIPTION: PETG SQUARE, MEDIA (LAB QUALITY) BOTTLES			