


| CATALOG NO. | NOMINAL CAPACITY<br>(oz) |   |    |    | MATERIALS |                |
|-------------|--------------------------|---|----|----|-----------|----------------|
|             |                          |   |    |    | BOTTLE    | CLOSURE        |
| 002100      | 4                        | 8 | 16 | 32 | FEP       | TEFZEL<br>ETFE |

| SIZE CODE.....:                      | -0004     | -0008      | -0016     |            |           |            |
|--------------------------------------|-----------|------------|-----------|------------|-----------|------------|
| SIZE (NOM VOL mL/oz)                 | 125/4     | 250/8      | 500/16    |            |           |            |
| NOM. BRIM CAP. (mL):                 | 145       | 260        | 520       |            |           |            |
| NECK FINISH (mm)....:                | 33-415    | 43-415     | 48-415    |            |           |            |
| THREAD (T.P.I.).....:                | 6         | 6          | 6         |            |           |            |
| NOM. WEIGHT<br>W/CLOSURE (gms).....: | 42        | 68         | 105       |            |           |            |
| A (NECK I.D. @ TOP):                 | 0.97 ±.04 | 24.6 ±1.0  | 1.28 ±.04 | 32.5 ±1.0  | 1.49 ±.05 | 37.8 ±1.3  |
| B (HT. W/CLOSURE)....:               | 4.61 ±.08 | 117.1 ±2.0 | 5.03 ±.08 | 127.8 ±2.0 | 6.55 ±.10 | 166.4 ±2.5 |
| C (HT. W/O CLOSURE):                 | 4.45 ±.08 | 113.0 ±2.0 | 4.86 ±.08 | 123.4 ±2.0 | 6.29 ±.10 | 159.8 ±2.5 |
| D (DIA).....:                        | 1.81 ±.06 | 46.0 ±1.5  | 2.34 ±.07 | 59.4 ±1.8  | 2.86 ±.08 | 72.6 ±2.0  |
| W.T.(NOM. WALL THK):                 | 0.04      | 1.0        | 0.04      | 1.0        | 0.04      | 1.0        |

| SIZE CODE.....:                      | -0032     |            |  |
|--------------------------------------|-----------|------------|--|
| SIZE (NOM VOL mL/oz)                 | 1000/32   |            |  |
| NOM. BRIM CAP. (mL):                 | 1060      |            |  |
| NECK FINISH (mm)....:                | 53-415    |            |  |
| THREAD (T.P.I.).....:                | 6         |            |  |
| NOM. WEIGHT<br>W/CLOSURE (gms).....: | 180       |            |  |
| A (NECK I.D. @ TOP):                 | 1.64 ±.06 | 41.7 ±1.5  |  |
| B (HT. W/CLOSURE)....:               | 8.1       | 206        |  |
| C (HT. W/O CLOSURE):                 | 7.81 ±.12 | 198.4 ±3.0 |  |
| D (DIA).....:                        | 3.54 ±.10 | 89.9 ±2.5  |  |
| W.T.(NOM. WALL THK):                 | 0.04      | 1.0        |  |

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



SPECS ON BACK

NOT TO SCALE

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

**ThermoFisher**  
SCIENTIFIC

Thermo Fisher Scientific  
75 Panorama Creek Drive  
Rochester, NY 14625-2385, USA

|                                    |     |       |         |                    |      |
|------------------------------------|-----|-------|---------|--------------------|------|
| MADE                               | PB  | APP'D | PWC     | PSS 00-2100-XXXX-1 | REV. |
| CK'D                               | KOW |       | 8/14/14 |                    | 0814 |
| DESCRIPTION:<br>WIDE-MOUTH BOTTLES |     |       |         |                    |      |