


CATALOG NO.	MATERIALS	
	GASKET	LINER
002101	SILICONE	TEFLON ( TFE )

	INCH	MM
A ( DIAMETER ) . . . . . :	4.00 ±.03	101.6 ±0.8
B ( THICKNESS ) . . . . . :	0.10 ±.02	2.5 ±0.5



SPECS ON BACK

NOT TO SCALE

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

**ThermoFisher**  
SCIENTIFIC

Thermo Fisher Scientific  
75 Panorama Creek Drive  
Rochester, NY 14625-2385, USA

MADE	JGJ	APP'D	PVC	PSS 00-2101-2200-1S	REV.	-
CK'D	PVC		4/24/96			
DESCRIPTION:				GASKET FOR WIDE MOUTH EP TOX BOTTLE CLOSURE		