


CATALOG NO.....	2114-0006		2110-0006	
BOTTLE MATERIAL:	HDPE		PPCO	
CLOSURE MATL....	PP		PP	
SIZE (NOM VOL mL/oz)	175/ 6		175/ 6	
NOM. BRIM CAP. (mL):	185		185	
NECK FINISH (mm):	38-415		38-415	
THREAD (T.P.I.):	6		6	
NOM. WEIGHT W/CLOSURE (gms):	33		33	
	INCH	MM	INCH	MM
A (NECK I.D.):	1.120 ±.025	28.4±0.6	1.120 ±.025	28.4±0.6
B (HT. W/CLOSURE):	4.17 ±.08	105.9±2.0	4.17 ±.08	105.9±2.0
C (HT. W/O CLOSURE):	4.04 ±.08	102.6±2.0	4.04 ±.08	102.6±2.0
D (SQUARE):	2.03 ±.05	51.6±1.3	2.13 ±.05	54.1±1.3
W.T.(NOM. WALL THK):	0.060	1.52	0.060	1.52

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



SPECS ON BACK

NOT TO SCALE

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

**ThermoFisher**  
SCIENTIFIC

Thermo Fisher Scientific  
75 Panorama Creek Drive  
Rochester, NY 14625-2385, USA

MADE	TL	APP'D	PSS 00-211X-0006	REV.
CK'D	PB	GRP 2/29/12		0212

DESCRIPTION:  
SQUARE, WIDE-MOUTH  
LAB QUALITY BOTTLES