


CATALOG NUMBER.....	2135-3803	2135-5303	2135-8303			
MATERIAL	CLOSURE.....	PP				
	INSERT.....	SILICONE - PLATINUM CURED				
	TUBING.....	SILICONE - PLATINUM CURED				
CLOSURE SIZE (mm).....	38-430	53B	83B			
TUBING SIZES (COMBINATION).....	(2) 1/8 I.D. (1) 1/4 I.D.	(2) 1/4 I.D. (1) 1/8 I.D.	(2) 1/4 I.D. (1) 3/8 I.D.			
	INCH	MM	INCH	MM	INCH	MM
A DIA.....	1.66±.03	42.2±0.8	2.58±.02	65.5±0.5	4.00±.03	101.6±0.8
B (CLOSURE HT).....	1.17±.03	29.7±0.8	1.65±.02	41.9±0.5	2.82±.03	71.6±0.8
H (INSERT HT).....	0.42±.03	10.7±0.8	0.67±.03	17.0±0.8	0.91±.04	23.1±1.0
TUBING LENGTHS (ABOVE INSERT).....	12	305	12	305	12	305
TUBING LENGTHS (BELOW INSERT)	1/8 I.D.	16	406	6	152	—
	1/4 I.D.	6	152	20	508	34
	3/8 I.D.	—	—	—	—	6
NOM WT. (gms).....	80	140	345			

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

**ThermoFisher**  
**SCIENTIFIC**

Thermo Fisher Scientific  
75 Panorama Creek Drive  
Rochester, NY 14625-2385, USA

MADE	MC	APP'D	PVC	PSS 00-2135-XX03	REV.	0501
CK'D	PVC		5/15/01			
DESCRIPTION: FLEXIBLE TOP WORKS™ SYSTEM SCREW CLOSURES WITH SILICONE INSERTS AND 3 PORTS						