

CATALOG NO.	NOMINAL CAPACITY (GAL.)			MATERIALS		
				BOTTLE	CLOSURE	SEAL RING
002202	1/2	1	2	LDPE	PP	—
002203	1/2	1	2	PPCO	PP	TPE***
002204	1/2	1	2	OPAQUE AMBER PPCO	OPAQUE AMBER PP **	TPE***
DS2205	1/2	1	2	PC	PP	TPE***


**NOTES:**

1. SIZE CODES -0210 AND -0250 ARE FOR CATALOG NUMBER DS2205.
- \* 2. 1 AND 2 GALLON SIZES HAVE MOLDED-IN SHOULDER LOOP (FOR ATTACHING IDENTIFICATION TAG) AS SHOWN.
- \*\* 3. THE CLOSURE ON CATALOG NUMBER 002204-0020 IS WHITE PP.
- \*\*\* 4. SEAL RING APPLIES ONLY TO -0020 SIZE FOR THESE CATALOGS.

SIZE CODE.....:	-0005/-0210		-0250		-0010*		-0020*	
SIZE (NOM VOL L/GAL):	2/0.5		2.5/0.5		4/1		8/2	
NOM. BRIM CAP. (mL):	2100		2275		4300		8000	
NECK FINISH (mm)....:	38-430		38-430		38-430		53B	
THREAD (T.P.I.).....:	6		6		6		4	
NOM. WEIGHT W/CLOSURE (gms).....:	180		195		370		700	
	INCHES	MM	INCHES	MM	INCHES	MM	INCHES	MM
A (NECK I.D.).....:	1.00 ±.05	25.4 ±1.3	1.00 ±.05	25.4 ±1.3	1.00 ±.05	25.4 ±1.3	1.46 ±.10	37.1 ±2.5
B (HT. W/CLOSURE)...:	9.80 ±.12	248.9 ±3.0	11.60 ±.14	294.6 ±3.6	13.10 ±.15	332.7 ±3.8	16.1 ±.3	409 ±8
C (HT. W/O CLOSURE):	9.60 ±.12	243.8 ±3.0	11.40 ±.14	289.6 ±3.6	12.90 ±.15	327.7 ±3.8	15.9 ±.2	404 ±5
D (DIA).....:	4.75 ±.10	120.7 ±2.5	4.75 ±.10	120.7 ±2.5	6.00 ±.12	152.4 ±3.0	7.60 ±.13	193.0 ±3.3
W.T.(NOM WALL THK):	0.07	1.8	0.08	2.0	0.10	2.5	0.12	3.0

**NOTES:**

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



**CONFIDENTIAL**

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

**ThermoFisher**  
**SCIENTIFIC**

Thermo Fisher Scientific  
75 Panorama Creek Drive  
Rochester, NY 14625-2385, USA

MADE	PB	APP'D	GRP	PSS 00-220X-XXXX-1	REV. 0513
CK'D	KOW		5/14/13		
DESCRIPTION:				LARGE, NARROW-MOUTH ROUND BOTTLES	