


CATALOG NO.	NOMINAL CAPACITY (oz)					MATERIALS		
						BOTTLE	CLOSURE	TUBULATION
002402	4	8	16	24	32	LOPE	PP	LOPE
002408	-	8	16	-	-	RED LOPE	RED PP	RED LOPE

SIZE CODE.....:	-0125		-0250		-0500	
SIZE (NOM VOL mL/oz)	125/4		250/8		500/16	
NOM. BRIM CAP. (mL):	145		290		580	
NECK FINISH (mm)....:	24-415		38-415		38-415	
THREAD (T.P.I.).....:	8		6		6	
NOM. WEIGHT W/CLOSURE (gms).....:	25		45		58	
	INCH	MM	INCH	MM	INCH	MM
A (NECK I.D.).....:	0.69 ±0.03	17.5 ±0.8	1.16 ±0.03	29.5 ±0.8	1.16 ±0.03	29.5 ±0.8
B (HT. W/CLOSURE)....:	4.70 ±0.08	119.4 ±2.0	5.42 ±0.08	137.7 ±2.0	6.81 ±0.09	173.0 ±2.3
C (HT. W/O CLOSURE):	4.62 ±0.08	117.3 ±2.0	5.31 ±0.08	134.9 ±2.0	6.69 ±0.09	170.0 ±2.3
D (DIA).....:	1.90 ±0.06	48.3 ±1.5	2.46 ±0.07	62.5 ±1.8	3.02 ±0.08	76.7 ±2.0
W.T.(NOM. WALL THK):	0.030	0.8	0.030	0.8	0.040	1.0

SIZE CODE.....:	-0750		-1000	
SIZE (NOM VOL mL/oz)	750/24		1000/32	
NOM. BRIM CAP. (mL):	840		1140	
NECK FINISH (mm)....:	38-415		43-415	
THREAD (T.P.I.).....:	6		6	
NOM. WEIGHT W/CLOSURE (gms).....:	82		105	
	INCH	MM	INCH	MM
A (NECK I.D.).....:	1.16 ±0.03	29.5 ±0.8	1.30 ±0.03	33.0 ±0.8
B (HT. W/CLOSURE)....:	8.05 ±0.10	204.5 ±2.5	8.25 ±0.10	209.6 ±2.5
C (HT. W/O CLOSURE):	7.92 ±0.10	201.2 ±2.5	8.12 ±0.10	206.2 ±2.5
D (DIA).....:	3.34 ±0.08	84.8 ±2.0	3.80 ±0.08	96.5 ±2.0
W.T.(NOM. WALL THK):	0.045	1.1	0.050	1.3

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



SPECS ON BACK

NOT TO SCALE

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

**ThermoFisher
SCIENTIFIC**

Thermo Fisher Scientific
75 Panorama Creek Drive
Rochester, NY 14625-2385, USA

MADE	MCC	APP'D	PVC	PSS 00-240X-XXXX-2	REV. 0593
CK'D	PVC		5/11/93		

DESCRIPTION:
WIDE MOUTH
UNITARY WASH BOTTLES