


CATALOG NO.	NOMINAL CAPACITY (mL/oz)		MATERIALS	
			BOTTLE	CLOSURE
002500	200/6	205/6	PSF	PP
002505	200/6	205/6	PPCO	PP
002507	200/6	205/6	PETG	HDPE

SIZE CODE.....:	-0280		-0380	
NOM. CAP. TO NECK(mL/oz):	200/6		210/6	
NOM. BRIM CAP. (mL):	210		220	
NECK FINISH (mm):	28-415		38-415	
THREAD (T.P.I.):	6		6	
NOM. WEIGHT W/CLOSURE (gms):	33		40	
	INCHES	MM	INCHES	MM
A (NECK I.D.):	0.84 ± 0.02	21.3 ± 0.5	1.14 ± 0.03	29.0 ± 0.8
B (HT. W/CLOSURE):	6.02 ± 0.09	152.9 ± 2.3	6.02 ± 0.03	152.9 ± 0.8
C (HT. W/O CLOSURE):	5.90 ± 0.09	149.9 ± 2.3	5.90 ± 0.09	149.9 ± 2.3
D (SQUARE):	1.75 ± 0.06	44.5 ± 1.5	1.75 ± 0.06	44.5 ± 1.5
W.T.(NOM. WALL THK):	0.060	1.52	0.060	1.52

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



SPECS ON BACK

NOT TO SCALE

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

ThermoFisher
SCIENTIFIC

Thermo Fisher Scientific
75 Panorama Creek Drive
Rochester, NY 14625-2385, USA

MADE	MCC	APP'D	PVC	PSS 00-250X-XXXX-1	REV.	0595
CK'D	PVC		5/8/95			

DESCRIPTION:

DILUTION BOTTLES