


CATALOG NO.....:		2600-0012
SIZE (NOM VOL (L)).....:		12
NOM. CAP. (L).....:		15
MATERIALS	BOTTLE.....:	PC
	CLOSURES.....:	WHITE PP
	GASKET.....:	SILICONE
NECK FINISH (mm)	TOP.....:	100
	PORTS.....:	38-430
THREAD (T.P.I.)	TOP.....:	5
	PORTS.....:	6
NOM. WEIGHT W/CLOSURES (gms).....:		1340
	INCH	MM
B (HT. W/ TOP CLOSURE).....:	16.7 ±0.3	424 ±8
C (HT. W/O TOP CLOSURE).....:	16.5 ±0.3	419 ±8
D (DIA).....:	11.4 ±0.2	290 ±5
I (NECK I.D.).....:	3.48 ±0.10	88.4 ±2.5
W.T.(NOM. WALL THK).....:	0.075	1.9



NOT TO SCALE

SPECS ON BACK

NOT TO SCALE

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

ThermoFisher
SCIENTIFIC

Thermo Fisher Scientific
75 Panorama Creek Drive
Rochester, NY 14625-2385, USA

MADE	MCC	APP'D	PVC	PSS 00-2600-0012	REV.	-
CK'D	PVC		11/3/92			
DESCRIPTION: 12L CULTURE VESSEL W/ PORTS						