


CATALOG NO.....:	3105-0015		3105-0050	
MATERIAL.....:	PC		PC	
SIZE (NOM VOL mL):	15		50	
NOM. BRIM CAP. (mL):	16		55	
NOM. WEIGHT(gms):	7		15	
	INCH	MM	INCH	MM
A (TOP I.D.).....:	0.59 ±0.01	15.0 ±0.3	1.04 ±0.01	26.4 ±0.3
B (BODY O.D.).....:	0.67 ±0.01	17.0 ±0.3	1.13 ±0.01	28.7 ±0.3
C (HEIGHT).....:	4.73 ±0.02	120.1 ±0.5	5.28 ±0.02	134.1 ±0.5



SPECS ON BACK

NOT TO SCALE

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

**ThermoFisher**  
SCIENTIFIC

Thermo Fisher Scientific  
75 Panorama Creek Drive  
Rochester, NY 14625-2385, USA

MADE	MCC	APP'D	PVC	PSS 00-3105-XXXX-1	REV.	-
CK'D	PVC		3/24/92			

DESCRIPTION:

CONICAL CENTRIFUGE TUBES