


CATALOG NO.	NOMINAL CAPACITY, (mL)				MATERIALS	
					TUBE	CLOSURE
3119	10	28	30	50 ¹	PPCO	PP

SIZE CODE.....	-0010		-0028		-0030		-0050 ¹	
NOM. CAP. (mL).....	10		28		30		50	
CLOSURE SIZE (mm)...	13		20		20		24	
	INCH	MM	INCH	MM	INCH	MM	INCH	MM
A (NECK I.D.).....	0.32	8.1	0.54	13.7	0.54	13.7	0.64	16.3
B (HT. W/CLOSURE)...	3.17 ±.03	80.5 ±.8	3.70 ±.03	94.0 ±.8	3.68 ±.03	93.5 ±.8	4.17 ±.045	105.9 ±1.1
C (TUBE O.D.).....	0.635 ^{+0.000} / _{-.018}	16.13 ^{+0.00} / _{-.46}	1.016 ^{+0.000} / _{-.040}	25.81 ^{+0.00} / _{-1.02}	1.016 ^{+0.000} / _{-.040}	25.81 ^{+0.00} / _{-1.02}	1.135 ^{+0.000} / _{-.040}	28.83 ^{+0.00} / _{-1.02}

¹ NOTE: ACTUAL BRIM CAPACITY OF CATALOG NO. 3119-0050 IS 42mL.

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



SPECS ON BACK

NOT TO SCALE

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

ThermoFisher
SCIENTIFIC

Thermo Fisher Scientific
75 Panorama Creek Drive
Rochester, NY 14625-2385, USA

MADE TL	APP'D GRP	PSS 00-3119-XXXX-1	REV. 1009
CK'D DV	10/23/09		

DESCRIPTION:

OAK RIDGE CENTRIFUGE TUBES