


CATALOG NUMBER.....		3139-0010		3139-0016	
NOM. CAP. (mL).....		10		16	
MATERIALS	TUBE.....	PPCO		PPCO	
	GASKET.....	RED SILICONE		RED SILICONE	
	CLOSURE.....	WHITE PP		WHITE PP	
CLOSURE SIZE (mm).....		13		15	
NOM. WEIGHT W/ GASKET & CLOSURE.....		4		7	
		INCH	MM	INCH	MM
A (NECK I.D.).....		0.32	8.1	0.38	9.6
B (TUBE O.D.).....		$.635^{+0}_{-.02}$	$16.13^{+0}_{-.5}$	$.715^{+0}_{-.03}$	$18.16^{+0}_{-.8}$
C (HT. W/O CLOSURE).....		3.09	78.5	3.83	97.3
D (HT. W/ CLOSURE).....		3.26	82.8	4.19	106.4

CATALOG NUMBER.....		3139-0030		3139-0050	
NOM. CAP. (mL).....		30		50	
MATERIALS	TUBE.....	PPCO		PPCO	
	GASKET.....	RED SILICONE		RED SILICONE	
	CLOSURE.....	WHITE PP		WHITE PP	
CLOSURE SIZE (mm).....		20		24	
NOM. WEIGHT W/ GASKET & CLOSURE.....		11		18	
		INCH	MM	INCH	MM
A (NECK I.D.).....		0.54	13.7	0.64	16.3
B (TUBE O.D.).....		$1.016^{+0}_{-.04}$	$25.81^{+0}_{-1.0}$	$1.135^{+0}_{-.04}$	$28.83^{+0}_{-1.0}$
C (HT. W/O CLOSURE).....		3.60	91.4	4.06	103.1
D (HT. W/ CLOSURE).....		4.02	102.1	4.47	113.5

NOTES:

- ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
- THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



CAT. NO. 3139-0010

CAT. NO. 3139-0016
3139-0030
3139-0050
(TYP)

SPECS ON BACK

NOT TO SCALE

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

ThermoFisher
SCIENTIFIC

Thermo Fisher Scientific
75 Panorama Creek Drive
Rochester, NY 14625-2385, USA

MADE	TL	APP'D	PSS 00-3139-00XX-1	REV.
CK'D	OV	GRP		1009
10/23/09				
DESCRIPTION: OAK RIDGE CENTRIFUGE TUBES WITH SEALING CAP				