


CATALOG NO.	NOMINAL CAPACITY, (mL)	MATERIALS	
		TUBE	CLOSURE
3148	35	PPCO	PP

SIZE CODE.....:	-0050	
NOM. CAP. (mL).....:	35	
CLOSURE SIZE (mm)....:	24	
	INCH	MM
A (NECK I.D.).....:	0.640	16.26
C (HT. W/CLOSURE)....:	4.09 ±.03	103.9 ±.8
D (TUBE O.D.).....:	1.115 ±.025	28.32 ±.64



SPECS ON BACK

NOT TO SCALE

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

**ThermoFisher**  
SCIENTIFIC

Thermo Fisher Scientific  
75 Panorama Creek Drive  
Rochester, NY 14625-2385, USA

MADE	TL	APP'D	GRP
CK'D	DV	10/23/09	

PSS 00-3148-0050

REV.  
1009

DESCRIPTION:  
CONICAL OAK RIDGE CENTRIFUGE TUBE