


CATALOG NO.....:	3650-0025	3650-0050	3650-0100			
NOM. CAP. (mL).....:	25	50	100			
WEIGHT (gms).....:	40	58	76			
CALIBRATION INCREMENTS (ml)....:	.01	.01	.02			
	INCH	MM	INCH	MM	INCH	MM
A (BODY I.D.).....:	.37	9.4	.42	10.7	.59	15.0
B (BODY O.D.).....:	.50	12.7	.55	14.1	.73	18.5
H (OVERALL HEIGHT)....:	20.8	528	28.3	719	28.4	721

COMPONENT	MATERIAL
BURET BODY	ACRYLIC
STOPCOCK ASSEMBLY	—
STOPCOCK & TIP	PMP
PLUG	TFE
O-RINGS	VITON
HANDLE	PP

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



SPECS ON BACK

NOT TO SCALE

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

**ThermoFisher**  
SCIENTIFIC

Thermo Fisher Scientific  
75 Panorama Creek Drive  
Rochester, NY 14625-2385, USA

MADE	MC <sup>2</sup>	APP'D	PVC	PSS 00-3650-0XXX-1	REV.	1100
CK'D	PVC		11/20/00			

DESCRIPTION:  
TRANSPARENT BURETS