

| CATALOG NO. | NOMINAL CAPACITY (oz) | | | | | | MATERIALS | |
|-------------|-----------------------|---|---|---|----|----|-----------|---------|
| | | | | | | | CYLINDER | *BASE |
| 003673 | 1 | 2 | 4 | 8 | 16 | 32 | PMP | BLUE PP |


* 1, & 2oz ARE ONE-PIECE CONSTRUCTION IN PMP;
4, 8, 16, & 32oz ARE TWO-PIECE CONSTRUCTION WITH A PP BASE.

| SIZE CODE.....: | -0001 | | -0002 | | -0004 | |
|------------------------|------------|------------|------------|------------|------------|------------|
| SIZE (NOM VOL mL/oz): | 30/1 | | 60/2 | | 125/4 | |
| BRIM CAP (mL).....: | 58 | | 113 | | 184 | |
| NOM. WEIGHT (gms)....: | 15 | | 26 | | 29 | |
| | INCH | MM | INCH | MM | INCH | MM |
| A (OD @ TOP).....: | 1.66 ±0.03 | 42.2 ±0.8 | 2.03 ±0.03 | 51.6 ±0.8 | 2.39 ±0.04 | 60.7 ±1.0 |
| B (DIM W/SPOUT).....: | 1.85 ±0.03 | 47.0 ±0.8 | 2.33 ±0.04 | 59.2 ±1.0 | 2.74 ±0.04 | 69.6 ±1.0 |
| C (ID @ TOP).....: | 1.58 ±0.03 | 40.1 ±0.8 | 1.94 ±0.03 | 49.3 ±0.8 | 1.95 ±0.03 | 49.5 ±0.8 |
| D | 4.20 ±0.06 | 106.7 ±1.5 | 5.37 ±0.08 | 136.4 ±2.0 | 5.84 ±0.08 | 148.3 ±2.0 |
| E (ACROSS FLATS)....: | 1.99 ±0.04 | 50.5 ±1.0 | 2.28 ±0.04 | 57.9 ±1.0 | 2.25 ±0.04 | 57.2 ±1.0 |

| SIZE CODE.....: | -0008 | | -0016 | | -0032 | |
|------------------------|------------|------------|------------|------------|-------------|------------|
| SIZE (NOM VOL mL/oz): | 250/8 | | 500/16 | | 1000/32 | |
| BRIM CAP (mL).....: | 380 | | 725 | | 1470 | |
| NOM. WEIGHT (gms)....: | 49 | | 94 | | 180 | |
| | INCH | MM | INCH | MM | INCH | MM |
| A (OD @ TOP).....: | 3.14 ±0.04 | 79.8 ±1.0 | 3.88 ±0.05 | 98.6 ±1.3 | 4.82 ±0.06 | 122.4 ±1.5 |
| B (DIM W/SPOUT).....: | 3.54 ±0.04 | 89.9 ±1.0 | 4.34 ±0.06 | 110.2 ±1.5 | 5.43 ±0.08 | 137.9 ±2.0 |
| C (ID @ TOP).....: | 3.02 ±0.04 | 76.7 ±1.0 | 3.73 ±0.05 | 94.7 ±1.3 | 4.50 ±0.06 | 114.3 ±1.5 |
| D | 7.10 ±0.09 | 180.3 ±2.3 | 9.27 ±0.10 | 235.5 ±2.5 | 11.34 ±0.10 | 288.0 ±2.5 |
| E (ACROSS FLATS)....: | 2.90 ±0.04 | 73.7 ±1.0 | 4.00 ±0.06 | 101.6 ±1.5 | 4.88 ±0.06 | 124.0 ±1.5 |

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



SPECS ON BACK

NOT TO SCALE

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

ThermoFisher
SCIENTIFIC

Thermo Fisher Scientific
75 Ponorama Creek Drive
Rochester, NY 14625-2385, USA

| | | | | | |
|---|--------------------|---------|-----|--------------------|------|
| MADE IN | NS/MC ² | APP'D | PVC | PSS 00-3673-XXXX-1 | REV. |
| CR'D | PVC | 5/23/02 | | | 0502 |
| DESCRIPTION: DOUBLE-SCALE PHARMACEUTICAL GRADUATES | | | | | |