


CATALOG NO.	NOMINAL CAPACITY, (mL)					MATERIALS	
						FLASK	CLOSURE
004001	50	100	250	500	1000	PMP	PP

SIZE CODE.....:	-0050		-0100		-0250	
NOM. CAP. (mL).....:	50		100		250	
CLOSURE SIZE (mm)....:	13		20		20	
NOMINAL WEIGHT WITH CLOSURE (gms).....:	15		20		40	
	INCH	MM	INCH	MM	INCH	MM
A (NECK I.D. AT TOP):	0.3	7.6	0.5	13	0.5	13
B (NECK O.D.).....:	0.46	11.7	0.6	15	0.7	18
C (BULB DIA).....:	1.82	46.2	2.18	55.4	3.0	76
D (HT W/O CLOSURE)....:	5.7	145	7.4	188	8.8	224
H (HT WITH CLOSURE)..:	5.8	147	7.5	191	8.9	226

SIZE CODE.....:	-0500		-1000	
NOM. CAP. (mL).....:	500		1000	
CLOSURE SIZE (mm)....:	24		28	
NOMINAL WEIGHT WITH CLOSURE (gms).....:	70		110	
	INCH	MM	INCH	MM
A (NECK I.D. AT TOP):	0.7	18	0.8	20
B (NECK O.D.).....:	0.8	20	0.9	23
C (BULB DIA).....:	3.7	94	4.7	119
D (HT W/O CLOSURE)....:	10.8	273.1	12.3	312
H (HT WITH CLOSURE)..:	10.9	269	12.4	315

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



SPECS ON BACK

NOT TO SCALE

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

**ThermoFisher**  
SCIENTIFIC

Thermo Fisher Scientific  
75 Ponorama Creek Drive  
Rochester, NY 14625-2385, USA

MADE	MC <sup>2</sup>	APP'D	PVC	PSS 00-4001-XXXX	REV.	0103
CR'D	PVC		1/13/03			

DESCRIPTION:

VOLUMETRIC FLASKS