



| CATALOG NUMBER | NOMINAL CAPACITY (mL) |     |     |      |      |      | MATERIALS |             |                   |
|----------------|-----------------------|-----|-----|------|------|------|-----------|-------------|-------------------|
|                |                       |     |     |      |      |      | BOTTLE    | CLOSURE     | LINER (70mm ONLY) |
| 004113         | 125                   | 250 | 500 | 1000 | 2000 | 2800 | PETG      | HDPE(WHITE) | LDPE/FOAM         |

| SIZE CODE.....:                    | -0125  |       | -0250  |       | -0500  |       |
|------------------------------------|--------|-------|--------|-------|--------|-------|
| NOM. CAP.....:                     | 125    |       | 250    |       | 500    |       |
| CLOSURE SIZE (mm).....:            | 38/430 |       | 38/430 |       | 45/430 |       |
| NOM. WEIGHT W/ CLOSURE (gms).....: | 35     |       | 45     |       | 94     |       |
|                                    | INCH   | MM    | INCH   | MM    | INCH   | MM    |
| A (NECK I.D.).....:                | 1.04   | 26.4  | 1.04   | 26.4  | 1.48   | 37.6  |
| B (HT. W/ CLOSURE).....:           | 4.44   | 112.8 | 5.39   | 136.9 | 7.00   | 177.8 |
| C (HT W/O CLOSURE):                | 4.34   | 110.2 | 5.24   | 133.1 | 6.83   | 173.5 |
| D (DIA.).....:                     | 2.61   | 66.3  | 3.25   | 82.6  | 4.00   | 101.6 |
| W.T. (NOM. WALL THK.):             | .03    | 0.7   | .03    | 0.7   | .05    | 1.3   |

| SIZE CODE.....:                    | -1000  |       | -2000  |       | -2800 |       |
|------------------------------------|--------|-------|--------|-------|-------|-------|
| NOM. CAP.....:                     | 1000   |       | 2000   |       | 2800  |       |
| CLOSURE SIZE (mm).....:            | 45/430 |       | 45/430 |       | 70    |       |
| NOM. WEIGHT W/ CLOSURE (gms).....: | 100    |       | 150    |       | 300   |       |
|                                    | INCH   | MM    | INCH   | MM    | INCH  | MM    |
| A (NECK I.D.).....:                | 1.40   | 35.6  | 1.40   | 35.6  | 2.40  | 61.0  |
| B (HT. W/ CLOSURE).....:           | 8.80   | 223.5 | 10.79  | 274.1 | 9.40  | 238.8 |
| C (HT W/O CLOSURE):                | 8.61   | 218.7 | 10.64  | 270.3 | 9.02  | 229.1 |
| D (DIA.).....:                     | 5.11   | 129.8 | 6.36   | 161.5 | 8.00  | 203.2 |
| W.T. (NOM. WALL THK.):             | .03    | 0.7   | .03    | 0.7   | .03   | 0.7   |

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



BOTTOM VIEW

SPECS ON BACK

NOT TO SCALE

**CONFIDENTIAL**

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

**ThermoFisher**  
SCIENTIFIC

Thermo Fisher Scientific  
75 Panorama Creek Drive  
Rochester, NY 14625-2385, USA

|      |     |        |                  |      |
|------|-----|--------|------------------|------|
| MADE | TL  | APP'D  | GRP              | REV. |
| CK'D | CHL | 3/5/13 | PSS 00-4113-XXXX | 0313 |

DESCRIPTION: **BAFFLED STERILE DISPOSABLE FLASKS**