


CAT. NO.	APPROX. SIZES (TOP I.D. IN MM)						MATERIAL
4256-	65	84	103	122	155	205	SEE CHART BELOW
SIZE CODE.....:	-0234		-0314		-0414		
MATERIAL.....:	PP		PP		PP		
NOM WT (gms).....:	10		18		29		
	INCH	MM	INCH	MM	INCH	MM	
A DIA.....:	2.69 ±0.04	68.3 ±1.0	3.44 ±0.04	87.4 ±1.0	4.30 ±0.06	109.2 ±1.5	
B (TOP I.D.).....:	2.45 ±0.04	62.2 ±1.0	3.15 ±0.04	80.0 ±1.0	3.87 ±0.04	98.3 ±1.0	
C DIA.....:	0.53 ±0.02	13.5 ±0.5	0.50 ±0.02	12.7 ±0.5	0.70 ±0.02	17.8 ±0.5	
D DIA.....:	0.39 ±0.02	9.9 ±0.5	0.37 ±0.02	9.4 ±0.5	0.44 ±0.02	11.2 ±0.5	
E (STEM HT).....:	0.82 ±0.02	20.8 ±0.5	0.85 ±0.02	21.6 ±0.5	1.47 ±0.04	37.3 ±1.0	
H (OVERALL HT)...	2.58 ±0.04	65.5 ±1.0	3.71 ±0.04	94.2 ±1.0	4.61 ±0.06	117.1 ±1.5	
SIZE CODE.....:	-0500		-0638		-0800		
MATERIAL.....:	HDPE		HDPE		HDPE		
NOM WT (gms).....:	41		73		95		
	INCH	MM	INCH	MM	INCH	MM	
A DIA.....:	4.99 ±0.06	126.7 ±1.5	6.3 ±0.1	160 ±3	8.0 ±0.1	203 ±3	
B (TOP I.D.).....:	4.62 ±0.06	117.3 ±1.5	5.9 ±0.1	150 ±3	7.6 ±0.1	193 ±3	
C DIA.....:	0.74 ±0.02	18.8 ±0.5	1.88 ±0.04	47.8 ±1.0	0.94 ±0.02	23.9 ±0.5	
D DIA.....:	0.51 ±0.02	13.0 ±0.5	1.03 ±0.02	26.2 ±0.5	0.78 ±0.02	19.8 ±0.5	
E (STEM HT).....:	1.81 ±0.04	46.0 ±1.0	2.61 ±0.04	66.3 ±1.0	1.79 ±0.04	45.5 ±1.0	
H (OVERALL HT)...	5.7 ±0.1	145 ±3	6.9 ±0.1	175 ±3	7.9 ±0.1	201 ±3	

- NOTES:
- ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
 - THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



NOT TO SCALE

SPECS ON BACK

NOT TO SCALE

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

ThermoFisher
SCIENTIFIC

Thermo Fisher Scientific
75 Panorama Creek Drive
Rochester, NY 14625-2385, USA

MADE	MCC	APP'D	PVC	PSS 00-4256-XXXX-1	REV.	1294
CK'D	PVC	12/28/94				

DESCRIPTION:

UTILITY FUNNELS