


CAT. NO.	APPROX. SIZES (FUNNEL TOP I.D. IN MM)					MATERIAL	
						UPPER	LOWER
4280-	59	75	99	132	175	PP	PP

SIZE CODE.....:	-0550		-0700		-0900	
NOM WT (gns).....:	23		55		80	
	INCH	MM	INCH	MM	INCH	MM
A DIA.....:	2.46 ±0.04	62.5 ±1.0	3.16 ±0.04	80.3 ±1.0	4.08 ±0.06	103.6 ±1.5
B (TOP I.D.).....:	2.30 ±0.04	58.4 ±1.0	2.96 ±0.04	75.2 ±1.0	3.88 ±0.06	98.6 ±1.5
C DIA.....:	0.55 ±0.02	14.0 ±0.5	0.66 ±0.02	16.8 ±0.5	0.74 ±0.02	18.8 ±0.5
D DIA.....:	0.41 ±0.02	10.4 ±0.5	0.44 ±0.02	11.2 ±0.5	0.48 ±0.02	12.2 ±0.5
E (STEM HT).....:	1.78 ±0.04	45.2 ±1.0	2.55 ±0.04	64.8 ±1.0	2.76 ±0.04	70.1 ±1.0
H (OVERALL HT)...	4.08 ±0.06	103.6 ±1.5	5.8 ±0.1	147.3 ±3	6.6 ±0.1	168 ±3

SIZE CODE.....:	-1100		-1500	
NOM WT (gns).....:	145		290	
	INCH	MM	INCH	MM
A DIA.....:	5.3 ±0.1	135 ±3	7.1 ±0.1	180 ±3
B (TOP I.D.).....:	5.2 ±0.1	132 ±3	6.9 ±0.1	175 ±3
C DIA.....:	1.14 ±0.04	29.0 ±1.0	1.15 ±0.04	29.2 ±1.0
D DIA.....:	0.87 ±0.02	22.1 ±0.5	0.80 ±0.02	20.3 ±0.5
E (STEM HT).....:	3.62 ±0.06	91.9 ±1.5	4.21 ±0.06	106.9 ±1.5
H (OVERALL HT)...	8.4 ±0.2	213 ±5	11.2 ±0.3	284 ±8



SPECS ON BACK

NOT TO SCALE

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

ThermoFisher
SCIENTIFIC

Thermo Fisher Scientific
75 Panorama Creek Drive
Rochester, NY 14625-2385, USA

MADE	MCC	APP'D	PVC	PSS 00-4280-XXXX-1	REV. 0195
CHK'D	PVC		1/18/95		

DESCRIPTION:

BÜCHNER FUNNELS