





CATALOG NUMBER.....:	5929-0016		5929-0020		5929-0030	
TEST TUBE SIZE.....:	16		20		30	
MATERIAL.....:	PC		PC		PC	
NOM. WEIGHT (gms)....:	250		230		260	
	INCH	MM	INCH	MM	INCH	MM
A (LENGTH).....:	11.28 ± .07	286.5 ± 1.8	11.28 ± .07	286.5 ± 1.8	12.8 ± .1	325.1 ± 2.5
B (WIDTH).....:	4.59 ± .05	116.6 ± 1.3	4.59 ± .05	116.6 ± 1.3	4.59 ± .05	116.6 ± 1.3
C (HEIGHT).....:	3.83 ± .05	97.3 ± 1.3	3.83 ± .05	97.3 ± 1.3	3.83 ± .05	97.3 ± 1.3
D (HEIGHT).....:	2.66 ± .03	67.6 ± 0.8	2.66 ± .03	67.6 ± 0.8	2.66 ± .03	67.6 ± 0.8
E (DIA.).....:	0.68 ± .01	17.3 ± 0.3	0.82 ± .01	20.8 ± 0.3	1.19 ± .01	30.2 ± 0.3
F.....:	1.00 ± .01	25.4 ± 0.3	1.00 ± .01	25.4 ± 0.3	1.44 ± .01	36.6 ± 0.3
G.....:	1.00 ± .01	25.4 ± 0.3	1.00 ± .01	25.4 ± 0.3	1.44 ± .01	36.6 ± 0.3
HOLE ARRAY.....:	4 x 10		4 x 10		3 x 8	



SPECS ON BACK

NOT TO SCALE

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

ThermoFisher
SCIENTIFIC

Thermo Fisher Scientific
75 Panorama Creek Drive
Rochester, NY 14625-2385, USA

MADE	PB	APP'D	GRP
CK'D	KOW	10/17/12	

PSS 00-5929-00XX

REV. 1012

TEST TUBE RACK