


CATALOG NO.....:	5976-0X13		5976-0X16	
TEST TUBE SIZE (mm):	13		16	
MATERIAL.....:	PP		PP	
NOM. WEIGHT (gms)...:	96		164	
	INCH	MM	INCH	MM
A (LENGTH).....:	7.9	201	9.8	249
B (WIDTH).....:	4.06	103.1	4.98	126.5
C (HEIGHT).....:	2.30	58.4	2.84	72.1
D (SQUARE).....:	0.54	13.7	0.69	17.5
HOLE ARRAY.....:	6 × 12		6 × 12	

COLOR	CATALOG NUMBER	
WHITE	5976-0013	5976-0016
BLUE	5976-0313	5976-0316
GREEN	5976-0413	5976-0416
RED	5976-0513	5976-0516

- NOTES:
- ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
 - THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



SPECS ON BACK NOT TO SCALE

CONFIDENTIAL
THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

ThermoFisher
SCIENTIFIC

Thermo Fisher Scientific
75 Panorama Creek Drive
Rochester, NY 14625-2385, USA

MADE IN CHINA	MC ² GRP	APP'D PVC 2/26/04	PSS 00-5976-0XXX	REV. 0204
DESCRIPTION: UNWIRE TEST TUBE RACKS				