


CATALOG NUMBER.....:	6379-0004	
NOM. CONTAINER CAPACITY (L).....:	4	
CLOSURE SIZE (mm).....:	38-430	
NOM. WT. (gms).....:	570	
	INCH	MM
FUNNEL TOP I.D.....:	5.5	140
NOM. OVERALL HT. W/ LID CLOSED.....:	17.6	447
NOM. OVERALL WIDTH W/ LID CLOSED.....:	7.0	178

COMPONENT	MATERIAL
COVER	RED HDPE
FUNNEL	HDPE
NECK SCREEN INSERT	PP
38-430 CLOSURE	PP
FUNNEL ADAPTER	PP
FILTER HOUSING/MEMBRANE	PP/PTFE
4L BOTTLE W/ CLOSURE	HDPE/PP



SPECS ON BACK

NOT TO SCALE

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

ThermoFisher
SCIENTIFIC

Thermo Fisher Scientific
75 Panorama Creek Drive
Rochester, NY 14625-2385, USA

MADE	MCC	APP'D	PVC	PSS 00-6379-0004	REV. 1101
CK'D	PVC	11/30/01			
DESCRIPTION: 4L SAFETY WASTE FUNNEL SYSTEM					