


CATALOG NUMBER.....:	6379-0010	
NOM. CONTAINER CAPACITY (L).....:	10	
CLOSURE SIZE (mm).....:	83B	
NOM. WT. (gms).....:	1920	
	INCH	MM
FUNNEL TOP I.D.....:	10	254
NOM. OVERALL HT. W/ LID CLOSED.....:	20.5	521
NOM. OVERALL WIDTH W/ LID CLOSED.....:	11.4	290

COMPONENT	MATERIAL
COVER	RED HDPE
FUNNEL	HDPE
NECK SCREEN INSERT	PP
83B CLOSURE	WHITE PP
FUNNEL ADAPTER	PP
FILTER HOUSING/MEMBRANE	PP/PTFE
10L CARBOY W/ CLOSURE	FLPE/WHITE PP



SPECS ON BACK

NOT TO SCALE

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

**ThermoFisher
SCIENTIFIC**

Thermo Fisher Scientific
75 Panorama Creek Drive
Rochester, NY 14625-2385, USA

MADE	MCC	APP'D	PVC	PSS 00-6379-0010	REV.	1101
CK'D	PVC	11/30/01				

DESCRIPTION: 10L SAFETY WASTE FUNNEL SYSTEM