


CATALOG NO.	NOMINAL CAPACITY (QTS)			MATERIALS	
				PAIL	HANDLE
007012	8	11	14	WHITE PP	BLACK NYLON-COATED WIRE

SIZE CODE.....:	-0080		-0110		-0140	
NOM. CAP. (L/QTS).....:	7.6/8		10.4/11		13.2/14	
NOM. BRIM CAP. (L).....:	8.4		12		15.2	
NOM. WEIGHT (gms).....:	430		550		680	
	INCH	MM	INCH	MM	INCH	MM
A (I.D. AT TOP).....:	8.75	222.3	9.74	247.4	10.62	269.8
B (HEIGHT).....:	10.2	259	11.5	292	12.4	315
C (HT. W/ HANDLE UP).....:	15.7	400	17.5	445	18.8	478
D (BASE O.D.).....:	7.75	196.9	8.65	219.7	9.16	240.3
E (LIP O.D. AT TOP).....:	9.5	241	10.5	267	11.4	290



NOT TO SCALE

SPECS ON BACK

NOT TO SCALE

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

**ThermoFisher**  
SCIENTIFIC

Thermo Fisher Scientific  
75 Panorama Creek Drive  
Rochester, NY 14625-2385, USA

MADE	MCC	APP'D	PVC	PSS 00-7012-0XXX-1	REV. 1195
CK'D	PVC		11/3/95		

DESCRIPTION:

BUCKETS