


CATALOG NO.	NOMINAL CAPACITY (QTS)			MATERIALS		
				PAIL	COVER	HANDLE
007102	8	11	14	* WHITE LOPE	WHITE LOPE	BLACK NYLON-COATED WIRE
* COLOR OF -0140 SIZE PAIL IS YELLOW.						
SIZE CODE.....	-0080		-0110		-0140	
NOM. CAP. (L/QTS).....	7.6/8		10.4/11		13.2/14	
NOM. BRIM CAP. (L).....	8.4		12		15.2	
NOM. WEIGHT WITH COVER (gms).....	550		700		900	
	INCH	MM	INCH	MM	INCH	MM
A (HT. W/O COVER).....	10.1	257	11.5	292	12.4	315
B (HT. W/ COVER).....	10.5	267	12.0	305	12.9	328
C (HT. W/ HANDLE UP).....	15.5	394	17.4	442	18.9	480
D (BASE O.D.).....	6.7	170	7.0	178	7.4	188
E (I.D. AT TOP).....	8.6	218	9.6	244	10.5	267
F (LIP O.D.).....	9.3	236	10.4	264	11.3	287
G (COVER O.D.).....	10.0	254	10.9	277	11.7	297
H (DIM. AT HANDLE'S WIDEST POINT).....	11.4	290	12.3	312	13.6	345

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



SPECS ON BACK

NOT TO SCALE

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

**ThermoFisher**  
SCIENTIFIC

Thermo Fisher Scientific  
75 Panorama Creek Drive  
Rochester, NY 14625-2385, USA

MADE	MCC	APP'D	PVC	PSS 00-7102-XXXX-1	REV.	—
CK'D	PVC		5/26/92			
DESCRIPTION: AIR-TIGHT PAILS						