


CATALOG NO.....	1060-85	
SIZE NOM VOL. (mL).....	2000	
NOMINAL BRIM CAP. (mL).....	2670	
MATERIALS	BOTTLE.....	PETG
	CLOSURE.....	TEAL HDPE
	MEMBRANE.....	PTFE
NECK FINISH.....	48mm DOUBLE START THREAD	
THREAD (T.P.I.).....	3	
NOM. WEIGHT W/ CLOSURE (gms).....	230	
	INCH	CM
A (NECK I.D. AT TOP).....	1.47	3.7
B (HT. W/O CLOSURE).....	10.67	27.1
C (HT. W/ CLOSURE).....	10.83	27.5
D (DIA).....	4.81	12.2
W.T. (NOM. WALL THK).....	0.03	0.08

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



SPECS ON BACK

NOT TO SCALE

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

ThermoFisher
SCIENTIFIC

Thermo Fisher Scientific
75 Panorama Creek Drive
Rochester, NY 14625-2385, USA

MADE	TL	APP'D	BR	PSS 1060-85	REV.
CK'D	JJ		3/27/13		0313
DESCRIPTION:					
ROLLER BOTTLE, PETG, 1.2X WITH VENTED CLOSURE					