


CATALOG NO.....	1060-05/-20	
SIZE NOM VOL. (mL).....	2000	
NOMINAL BRIM CAP. (mL).....	2670	
MATERIALS	BOTTLE.....	PETG
	CLOSURE.....	WHITE HDPE
NECK FINISH.....	48mm DOUBLE START THREAD	
THREAD (T.P.I.).....	3	
NOM. WEIGHT W/ CLOSURE (gms).....	230	
	INCH	CM
A (NECK I.D. AT TOP).....	1.47	3.7
B (HT. W/O CLOSURE).....	10.67	27.1
C (HT. W/ CLOSURE).....	10.83	27.5
D (DIA).....	4.81	12.2
W.T. (NOM. WALL THK).....	0.03	0.08

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



SPECS ON BACK

1060-05/20

NOT TO SCALE

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

**ThermoFisher**  
SCIENTIFIC

Thermo Fisher Scientific  
75 Panorama Creek Drive  
Rochester, NY 14625-2385, USA

MADE	TL	APP'D	BR	PSS 1060-XX	REV. 0313
CK'D	JJ		3/27/13		

DESCRIPTION:

ROLLER BOTTLE, PETG, 1.2X