


CATALOG NUMBER	SIZE (GAL.)	Ø A	B ±1/2"	Ø C	D ±1/4"	Ø E	WALL THK.
11100-0005	5	12-9/32	15-1/8	11	1	12-25/32	3/16
11100-0007	7-1/2	13 13/16	18	12 7/16	1 1/4	14 1/4	3/16
11100-0010	10	14 3/4	20 3/8	13 3/8	1 1/4	15 3/16	3/16
11100-0015	15	15	27 5/8	13 5/8	1 1/4	15 7/16	3/16
11100-0030	30	20	29 7/8	18 3/8	1 1/4	20 7/16	3/16
11100-0055	55	23 5/8	36 1/8	22	1 1/4	24 3/8	1/4
11100-0080	80	25 9/16	48 1/16	23 15/16	2	26 1/4	1/4
11100-0100	100	29 7/8	44 1/4	28 3/8	2	30 3/8	1/4
11100-0150	150	33 3/16	49 5/16	31 5/8	2	34 1/16	1/4
11100-0200	200	38 5/16	51	36 11/16	2	39 1/16	1/4



SPECS ON BACK

NOT TO SCALE

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

ThermoFisher
SCIENTIFIC

Thermo Fisher Scientific
75 Panorama Creek Drive
Rochester, NY 14625-2385, USA

MADE	PB	APP'D	GRP	PSS 11100-0XXX	REV.	0813
CK'D	RCB	8/20/13				
DESCRIPTION: CYLINDRICAL TANK ASSEMBLIES - LLDPE						