



NUMBER OF WELLS.....		1	
MATERIAL	DISH.....	POLYSTYRENE	
	LID.....	POLYSTYRENE	
		INCH	MM
A (DISH BASE LENGTH).....		5.03±.02	127.8±.5
B (DISH BASE WIDTH).....		3.37±.02	85.5±.5
C (DISH HEIGHT).....		0.57±.02	14.5±.5
D (LID LENGTH).....		5.04±.02	128.0±.5
E (LID WIDTH).....		3.38±.02	85.9±.5
F (LID HEIGHT).....		0.36±.02	9.1±.5
G (OVERALL HEIGHT).....		0.68±.05	17.3±1.3
H (WELL BOTTOM LENGTH, FLAT).....		4.73±.02	121.4±.5
J (WELL BOTTOM WIDTH, FLAT).....		3.07±.02	78.0±.5
K (WELL DEPTH).....		0.45±.02	11.6±.5
L (WELL BOTTOM THICKNESS).....		0.04±.02	1.0±.5


CATALOG NUMBER	DISH COLOR	TREATED
140156	CLEAR	YES
140175	WHITE	YES
267060	CLEAR	NO




TOP VIEW




FRONT VIEW



WIDTH CROSS-SECTION



LENGTH CROSS-SECTION



CORNER DETAIL

SPECS ON BACK

NOT TO SCALE

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

ThermoFisher
SCIENTIFIC

Thermo Fisher Scientific
75 Panorama Creek Drive
Rochester, NY 14625-2385, USA

MADE BY PELYCHATY	APP'D PWC	PSS 140156-S	REV. 0915
CK'D BY J. CLEMONS	9/21/15		

1-WELL DISH W/ LID