


NUMBER OF WELLS.....:		1	
MATERIAL	DISH.....:	POLYSTYRENE	
	LID.....:	POLYSTYRENE	
		INCH	MM
A (DISH BASE LENGTH).....:		5.03	127.8
B (DISH BASE WIDTH).....:		3.37	85.5
C (DISH HEIGHT).....:		0.57	14.5
D (LID LENGTH).....:		5.04	128.0
E (LID WIDTH).....:		3.38	85.9
F (LID HEIGHT).....:		0.36	9.1
G (OVERALL HEIGHT).....:		0.68	17.3
H (WELL BOTTOM LENGTH, FLAT).....:		4.78	121.4
J (WELL BOTTOM WIDTH, FLAT).....:		3.07	78.0
K (WELL DEPTH).....:		0.45	11.6
L (WELL BOTTOM THICKNESS).....:		.040	1.0

CATALOG NUMBER	DISH COLOR	TREATED
140156	CLEAR	YES
140175	WHITE	YES
267060	CLEAR	NO



SPECS ON BACK

NOT TO SCALE

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

**ThermoFisher**  
SCIENTIFIC

Thermo Fisher Scientific  
75 Panorama Creek Drive  
Rochester, NY 14625-2385, USA

MADE BY PB BR	APP'D BR 4/19/10	PSS 140156	REV. 0410
DESCRIPTION: 1-WELL DISH W/ LID			