




CATALOG NUMBER.....:		156800	
RECOMMENDED WORKING VOL...:		150	
MAXIMUM CAPACITY IN UPRIGHT POSITION (mL).....:		800	
CAP TYPE.....:		FILTER	
CULTURE AREA (cm ²).....:		225	
MATERIAL	FLASK.....:	PS	
	CAP.....:	WHITE HDPE	
		INCH	MM
A (NECK I.D.).....:	1.04	26.4	
B (LENGTH W/O CAP).....:	9.06	230.1	
C (LENGTH W/ CAP).....:	9.13	231.9	
D (WIDTH).....:	5.54	140.7	
E (HEIGHT).....:	1.57	39.9	



TOP VIEW



SIDE VIEW

END VIEW

SPECS ON BACK

NOT TO SCALE

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

ThermoFisher
SCIENTIFIC

Thermo Fisher Scientific
75 Panorama Creek Drive
Rochester, NY 14625-2385, USA

MADE	PB	APP'D			
CK'D					

PSS 156800

REV. 1107

DESCRIPTION: NUNC FLASKS NON-TREATED - STERILE