



CATALOG NUMBER	SUGGESTED WORKING VOL. (μL)	NUMBER OF WELLS	MATERIALS	
			TRAY	LID
163118	8	60	PS	PS

	INCH	MM
A (LENGTH).....:	3.31	84.1
B (WIDTH).....:	2.29	5.82
C (LENGTH).....:	3.37	85.6
D (WIDTH).....:	2.35	59.7
E (HEIGHT).....:	0.44	11.2
G (LID HEIGHT).....:	0.13	3.3
F (DIA.).....:	0.15	3.8

- ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
- THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



SPECS ON BACK

NOT TO SCALE

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

ThermoFisher
SCIENTIFIC

Thermo Fisher Scientific
75 Panorama Creek Drive
Rochester, NY 14625-2385, USA

MADE	P.B.	APP'D	JMC	PSS 163118	REV. 0309
CK'D	JMC	3/17/09			
DESCRIPTION:					
NUNC MINITRAY NUNCLON, STERILE					