


## Technical specifications

### Cell Factory 10 (CF10)

Material	Polystyrene
Dimensions	333 x 204 x 210 mm
Culture area	6,320 cm <sup>2</sup>
Surface	Nunculon™ Δ treated

Sterilization	Irradiated
Appearance	Optically clear
Volume markings	From 100 to 300 ml in 50 ml grads
Weight	2200 g



Reference on drawin	Description	Dimensions (mm)
A	Top length	333
B	Bottom length	329
C	Top depth	204
D	Bottom depth	200
E	Overall height of trays	188
F	Overall height	210

For instruction of handling CF10 please see chapter "Instructions";

**NUNC™**