



CATALOG NO.	COLOR	UNITS PER (PK/CASE)
165218	CLEAR	10/60
242811	CLEAR	10/60
264728	CLEAR	10/90
465219	CLEAR	10/60

TOTAL VOLUME (mL)	90	
SUGGESTED WORKING VOL. (mL)	35	
NOMINAL CULTURE AREA (cm ²)	84	
MATERIAL (PLATE & LID)	POLYSTYRENE	
	INCH	MM
A (TRAY LENGTH)	5.030±.02	127.76±.5
B (TRAY WIDTH)	3.365±.02	85.47±.5
C (TRAY LENGTH)	4.87±.02	123.7±.5
D (TRAY WIDTH)	3.20±.02	81.3±.5
E (FLANGE HEIGHT)	0.10±.02	2.5±.5
F (TRAY HEIGHT)	0.57±.02	14.5±.5
G (WELL DEPTH)	0.46±.02	11.7±.5
H (LID LENGTH)	5.02±.02	127.5±.5
J (LID WIDTH)	3.36±.02	85.3±.5
K (LID HEIGHT)	0.31±.02	7.9±.5
L (SINGLE ASSEMBLY HEIGHT)	0.63±.05	16.0±1.3
M (STACKED ASSEMBLY HEIGHT)	1.26±.1	32.0±2.54



SPECS ON BACK

NOT TO SCALE

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

ThermoFisher
SCIENTIFIC

Thermo Fisher Scientific
75 Panorama Creek Drive
Rochester, NY 14625-2385, USA

MADE PMP APP'D PVC
CKD JC 10/21/15

PSS 165218

REV. 1015

DESCRIPTION:

OMNITRAY