


NUMBER OF WELLS.....		4	
MATERIAL	DISH.....	POLYSTYRENE	
	LID.....	POLYSTYRENE	
		INCH	MM
A (DISH BASE LENGTH).....		5.03	127.8
B (DISH BASE WIDTH).....		3.37	85.5
C (DISH HEIGHT).....		0.57	14.5
D (LID LENGTH).....		5.04	128.0
E (LID WIDTH).....		3.38	85.9
F (LID HEIGHT).....		0.36	9.1
G (OVERALL HEIGHT).....		0.68	17.3
H (WELL BOTTOM WIDTH, FLAT)...		1.10	27.9
J (WELL BOTTOM LENGTH, FLAT)...		3.07	78.0
K (WELL DEPTH).....		0.45	11.6
L (BOTTOM THICKNESS).....		0.045	1.14

CATALOG NUMBER	TREATED
167063	YES
267061	NO



SPECS ON BACK

NOT TO SCALE

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

**ThermoFisher**  
SCIENTIFIC

Thermo Fisher Scientific  
75 Panorama Creek Drive  
Rochester, NY 14625-2385, USA

MADE	MC3JDD	APP'D	PVC	PSS 167063	REV.	0706
CK'D	PVC		1/15/02			
DESCRIPTION: 4-WELL DISH w/ LID						