



NUMBER OF WELLS.....:		8	
MATERIAL	DISH.....:	POLYSTYRENE	
	LID.....:	POLYSTYRENE	
		INCH	MM
A (DISH BASE LENGTH).....:		5.03	127.8
B (DISH BASE WIDTH).....:		3.37	85.5
C (DISH HEIGHT).....:		0.57	14.5
D (LID LENGTH).....:		5.04	128.0
E (LID WIDTH).....:		3.38	85.9
F (LID HEIGHT).....:		0.36	9.1
G (OVERALL HEIGHT).....:		0.68	17.3
H (WELL BOTTOM WIDTH, FLAT)...:		1.10	27.9
J (WELL BOTTOM LENGTH, FLAT)..:		1.48	37.6
K (WELL DEPTH).....:		0.45	11.6
L (BOTTOM THICKNESS).....:		0.048	1.22


CATALOG NUMBER	TREATED
167064	YES
267062	NO




TOP VIEW




FRONT VIEW



WIDTH CROSS-SECTION



LENGTH CROSS-SECTION



CORNER DETAIL

SPECS ON BACK

NOT TO SCALE

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

ThermoFisher
SCIENTIFIC

Thermo Fisher Scientific
75 Panorama Creek Drive
Rochester, NY 14625-2385, USA

MADE	MC ² JDD	APP'D	PSS 167064	REV. 0706
CK'D	PVC	PVC 2/25/06		

DESCRIPTION: 8-WELL DISH w/ LID