


CATALOG NO.....	1760-20/-90	
SIZE NOM VOL. (mL).....	1000	
NOMINAL BRIM CAP. (mL).....	2450	
MATERIALS	BOTTLE.....	PETG
	CLOSURE.....	WHITE HDPE
NECK FINISH.....	48mm DOUBLE START THR'D	
THREAD (T.P.I.).....	3	
NOM. WEIGHT W/ CLOSURE (gms).....	235	
	INCH	CM
A (NECK I.D. AT TOP).....	1.47	3.7
B (HT. W/O CLOSURE).....	10.67	27.1
C (HT. W/ CLOSURE).....	10.83	27.5
D (DIA).....	4.81	12.2
W.T. (NOM. WALL THK).....	0.02	0.05

2. This is not a production drawing, but is for informational purposes only, and is subject to change without notice.

1. All dimensions are nominal and in inches and millimeters unless otherwise noted.



SPECS ON BACK

1760-05/20

NOT TO SCALE

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

ThermoFisher
SCIENTIFIC

Thermo Fisher Scientific
75 Panorama Creek Drive
Rochester, NY 14625-2385, USA

MADE TL	APP'D BR	PSS 1760-XX	REV. 0313
CK'D JJ	3/27/13		
DESCRIPTION: ROLLER BOTTLE, PETG, 2X			