


CATALOG NUMBER	NOMINAL CAPACITY (L)	MATERIALS		
		BOTTLE	CLOSURE	CLOSURE LINER
2101-2200	2.2	TEFLON FEP	WHITE PP	TEFLON

NOMINAL VOLUME (L)	2.2	
NOMINAL BRIM CAPACITY (mL)	2280	
NECK FINISH	100	
THREAD TPI	5	
NOMINAL WEIGHT W/CLOSURE (g)	395	
	INCH	MM (REF)
A (NECK ID)	3.38±.12	85.9
B (HEIGHT W/CLOSURE)	9.4±.2	239
C (HEIGHT W/O CLOSURE)	9.12±.15	231.6
D (DIAMETER)	4.77±.08	121.2
E (NOMINAL WALL THICKNESS)	.040	1.02
F (MIN WALL THICKNESS AT MIDPOINT)	.030	.76

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



SPECS ON BACK

NOT TO SCALE

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

ThermoFisher
SCIENTIFIC

Thermo Fisher Scientific
75 Panorama Creek Drive
Rochester, NY 14625-2385, USA

MADE P. PELYCHATY	APP'D 4/22/16	PSS 2101-2200	REV. 0416
CK'D PWC	PWC		
DESCRIPTION: WIDE MOUTH EP TOX BOTTLE			