


CATALOG NO.....	2160-05/-20	
SIZE NOM VOL. (mL).....	1000	
NOMINAL BRIM CAP. (mL).....	2500	
MATERIALS	BOTTLE.....	PETG
	CLOSURE.....	WHITE HDPE
NECK FINISH.....	48mm DOUBLE START THR'D	
THREAD (T.P.I.).....	3	
NOM. WEIGHT W/ CLOSURE (gms).....	235	
	INCH	CM
A (NECK I.D. AT TOP).....	1.47	3.7
B (HT. W/O CLOSURE).....	10.67	27.1
C (HT. W/ CLOSURE).....	10.83	27.5
D (DIA).....	4.81	12.2
W.T. (NOM. WALL THK).....	0.02	0.05



SPECS ON BACK

NOT TO SCALE

NOTES:

1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

ThermoFisher
SCIENTIFIC

Thermo Fisher Scientific
75 Panorama Creek Drive
Rochester, NY 14625-2385, USA

MADE	PB	APP'D	PVC	REV.
CR'D	PVC	11/21/01	PSS 2160-XX	1101

DESCRIPTION:
ROLLER BOTTLE, PETG, 2.5XL